

SL-OCT can image in any cross-section from horizontal to vertical. The slit beam (on the patient's eye and/or in the video image) shows the location of the OCT scan.

**For quality white-to-white images, it is important that:**

- Both angles are visible – *to adjust: move slit lamp base right/left*
- Cornea is aligned with the green targeting arc – *to adjust: move device toward/away from patient*
- The image is level – *to adjust: move inner slit-lamp arm slightly right/left*

### Preparation (Section 9.2)

Start system and open *Image Acquisition* window. Be sure patient is:


- Properly positioned
- Sitting in a stable position
- Able to see the fixation target (Section 9.3)

### Align Image (Section 9.4)

- **Looking at the patient's eye**, adjust the slit lamp beam sharp and narrow over the apex of the cornea. When correctly positioned the slit beam is visible as a bright reflex from the pupil.
- **Looking at the screen**, adjust image by moving slit lamp base.
- **With Video Kit** – move slit lamp base right/left and toward/away from patient to align the pairs of white lights around the green target circles (see iris image above).

### Fine Tune Image (Section 9.4)

- Twist the joystick SLOWLY to bring slit beam to apex of cornea. As the beam approaches the apex you will see *bright reflex spots* appearing on the anterior surface of the cornea (see sample image below).
- At this point, ask the patient not to blink or move. Continue fine-tuning with the joystick until the reflex appears on both outer and inner corneal surfaces and the sharp reflex beam is visible through the anterior chamber (see acquisition window image above, left).



### Save Images (Section 7.4)

- When the desired number of “good” (green checkmark) images is visible click [Stop].
- Click an image once to view it in the large window; double-click an image to mark for saving.
- Select and save the best 3-5 images, or re-image this eye. Good scans:
  - > Have the green checkmark
  - > Are free of motion artifact (wavy cornea)
  - > Have bright reflex on both corneal surfaces and through image (see sample image at left).
  - > Are well centered
- After saving, scan the other eye or exit the acquisition window

*For additional imaging guidance refer to your SL-OCT Operating Instructions and SL-OCT Video Kit Operations Manual.*