
VISULAS 532s

Gebrauchsanweisung
User Manual
Mode d'emploi
Instrucciones de uso

ZEISS

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Knowledge of this user manual is required for operation of the instrument. You should therefore familiarise yourself with the contents of this manual and pay special attention to instructions concerning the safe operation of the instrument.

We reserve the right to make changes to the product in light of technical developments; these operating instructions will not be updated.

- © Unless expressly authorised, dissemination or duplication of this document and commercial exploitation or communication of its contents is not permitted. Persons in contravention of this copyright are liable to pay compensation for damages.
All rights reserved in the event of granting of patents or -registration as a utility patent.

	Page
Copyright	1
Contents	3
Figures	6
Safety notes	8
General	8
Standards and regulations	9
Notes on installation and use	10
Laser safety	13
Instrument safety systems	17
Warning and information labels	18
Description	22
Intended use	22
Operating principle	23
VISULAS 532s system design	24
VISULAS YAG III Combi system design	25
Controls, display elements and connectors	26
Laser slit lamp	26
Instrument base with headrest	28
Beam delivery in the laser slit lamp	30
Beam delivery in the YAG laser head (for VISULAS YAG III Combi)	31
Laser console	32
Control panel	34
VISULINK 532/U optical system	35
VISULINK 532/U with accessories	36
LIO 532 optical system	37
LIO 532 laser indirect ophthalmoscope with accessories	38
Operation	40
Installation of the instrument	40
Moving the IT 4L instrument table with the instrument	42
Menu navigation	43
VISULAS 532s description of menus	44
Power-up mode	44
Diagnosis mode	44
Main menu	45
System information	47
Program definition mode	48
Treatment mode	49

	Page
Operation with laser slit lamp	51
Adjusting the slit lamp and eyepieces	51
Positioning the contact lens	51
Laser treatment with VISULAS 532s and LSL	52
Positioning and focusing the laser spot	52
Laser configuration	52
Operation with the VISULINK 532/U	53
Mounting the VISULINK 532/U on Carl Zeiss Meditec- and Haag-Streit slit lamps	53
Removing the VISULINK 532/U from the slit lamp	54
Adjusting the slit lamp and the VISULINK 532/U	55
Laser treatment using the VISULINK 532/U	57
Operation with the LIO 532 laser indirect ophthalmoscope ..	58
Connecting the cold light source	58
Adjusting the head strap	58
Adjusting the binocular body	59
Examining the fundus	59
Replacing the +2 D lenses	60
Laser treatment with the LIO 532	60
Preparation	60
Laser treatment	61
Operation with the intraocular probe (Endoprobe)	62
Description	62
Physician's safety filter	63
Application	64
Sterile operation	65
Laser treatment	66
Maintenance, Miscellaneous	68
Troubleshooting tables	68
Warning messages and notes	68
Faults during operation of the LIO 532	70
Faults during operation of the laser slit lamp	71
Faults during operation of the VISULINK 532/U	72
Replacing the halogen bulb in the slit lamp	73
Care of the instrument	74
Transportation instructions	76
Disposal	76
Disposal of the product within the EU	76
Safety inspections	77

	Page
Calibrating the power measurement system	78
Calibration procedure	79
Calibrating the laser console	79
Calibrating the applicators	80
Ordering information, accessories and spare parts	81
Technical specifications	83
VISULAS 532s laser system	83
LSL 532s laser slit lamp	84
LSL YAG III Combi laser slit lamp	84
VISULINK 532/U	84
LIO 532 laser indirect ophthalmoscope	84
Manufacturer's Declaration	85
Important safety information	Page
Safety notes	8
Safety inspections	77

Fig. 1	Warning and information labels on the laser slit lamp LSL 532s	18
Fig. 2	Warning and information labels on the laser slit lamp LSL YAG III Combi	19
Fig. 3	Warning and information labels on the laser console ..	20
Fig. 4	Warning and information labels on the LIO 532 and VISULINK 532/U	21
Fig. 5	Generation of laser radiation with $\lambda = 532 \text{ nm}$	23
Fig. 6	VISULAS 532s components	24
Fig. 7	VISULAS YAG III Combi components	25
Fig. 8	LSL 532s laser slit lamp	27
Fig. 9	LSL YAG III Combi laser slit lamp	27
Fig. 10	Instrument base with headrest	29
Fig. 11	Beam delivery in the LSL 532s laser slit lamp	30
Fig. 12	Beam delivery in the YAG laser head	31
Fig. 13	Laser console with control panel	32
Fig. 14	Connectors and controls on the laser console	33
Fig. 15	Control panel, guide rod swivelled out	34
Fig. 16	Removing the control panel	34
Fig. 17	VISULINK 532/U optical system	35
Fig. 18	VISULINK 532/U with accessories	36
Fig. 19	Operating principle	37
Fig. 20	Controls on the LIO 532 laser indirect ophthalmoscope and the cold light source	38
Fig. 21	IT 4L with VISULAS 532s	40
Fig. 22	IT 4L with VISULAS YAG III Combi	40
Fig. 23	General menu navigation	43
Fig. 24	Power-up mode	44
Fig. 25	Diagnosis mode	44
Fig. 26	Main menu	45
Fig. 27	System information	47
Fig. 28	Programme definition mode	48
Fig. 29	Treatment STANDBY mode; a) for the LSL and VISULINK, b) for the LIO and ENDO	49
Fig. 30	Treatment READY mode; a) for the LSL and VISULINK, b) for the LIO and ENDO	50
Fig. 31	Removing the VISULINK 532/U	54
Fig. 32	VISULINK 532/U mounted on Haag-Streit slit lamp and Carl Zeiss Meditec SL 130	54
Fig. 33	Adjusting the laser spot	56
Fig. 34	Mounting the focusing aid on the VISULINK 532/U ...	56
Fig. 35	Replacing the lenses	60
Fig. 36	Intraocular probe	62
Fig. 37	Intraocular probe	63
Fig. 38	Control panel with sterile stylus and holder	65
Fig. 39	Confirmation dialog box for physician's safety filter type	66
Fig. 40	Replacing the halogen bulb	72

General

This instrument has been developed and tested in accordance with national and international regulations. A high degree of instrument safety is thus ensured.

The following chapter contains a summary of the most important safety information.

All safety notes and information in this manual and on the labels on the instrument must be observed.

Warning

Danger to the user!

Caution

Danger of damage to the instrument!

Symbol on type label

Type B medical device in accordance with DIN EN 60601-1.

Disconnect the instrument from the power supply before servicing.

Warning - laser radiation!

Laser radiation emission at end of optical fibre!

Note

Information and notes for a better understanding of the operating procedures during operation of the instrument

Note for VISULAS YAG III Combi

Special note when using the VISULAS YAG III Combi

Warning

Correct operation of the instrument is imperative for its safe function. You should therefore ensure that you are thoroughly familiar with these operating instructions before using the instrument.

Laser damage may occur if the instrument is operated or used for the treatment of patients in a manner other than that described in these operating instructions.

Operating instructions for other equipment used with this instrument must be observed.

Further information may be obtained from Carl Zeiss Meditec authorised specialists or from authorised representatives.

Standards and regulations

- ❑ Carl Zeiss Meditec works in accordance with a certified quality management system.
- ❑ This instrument is equipped with a power display, a key switch, remote interlock contact and all required warning and information labels in accordance with the relevant standards.
- ❑ All relevant national accident prevention regulations for this instrument must be observed.
- ❑ In some countries (e.g. in the USA), national regulations specify that this instrument may only be used under the supervision of a medical doctor.

Warning

In the USA, this device may only be sold to a medical doctor or on the basis of an order from a medical doctor.

- ❑ This instrument is a class IV (4) laser device. The safety regulations applying to this laser classification must be observed.
- ❑ This instrument is a class IIb medical instrument as defined by the European Medical Device Directive (MDD).
- ❑ The instrument meets the requirements of the EU- Medical Device Directive 93/42/EEC and its national implementation in the form of the German Medical Product Act (--> *Manufacturer's Declaration on p.85*).
- ❑ To ensure proper function of the instrument, regular safety inspections should be carried out. The instrument should be inspected by a Carl Zeiss Meditec authorised specialist annually and the results entered in the instrument logbook. See also the chapter on *Safety inspections*, p.77 of these operating instructions.
- ❑ An instrument logbook must be kept.

Notes on installation and use

General

- ❑ This instrument may only be operated by appropriately qualified and trained personnel. It is the responsibility of the institution operating the instrument to train and instruct operating personnel.
- ❑ People working in the laser hazard area must be given instruction concerning safety regulations, safety precautions and operation of the instrument at least once a year. A written record of this instruction, listing the participants, must be kept.
- ❑ The operating instructions and instrument logbook should be readily accessible to persons operating the instrument at all times.
- ❑ The instrument should be used for the applications described only. Damage occurring through the instrument being used in a manner other than that intended is the responsibility of the user.
- ❑ The LIO 532 laser indirect ophthalmoscope may only be used as an applicator in combination with the Carl Zeiss Meditec VISULAS 532s laser console.
- ❑ The VISULINK 532/U may only be used as an applicator in combination with the Carl Zeiss Meditec VISULAS 532s laser console.
- ❑ The instruments supplied must not be operated
 - in areas in which there is a danger of explosion
 - or in the presence of volatile anaesthetics or flammable solvents, such as alcohol, benzene, etc.
- ❑ Do not store or use this device in damp rooms. Do not allow water to drip or splash near the instrument.
- ❑ Do not place fluid-filled containers on top of the laser console.
- ❑ Switch off the instrument immediately if smoke, sparks or unusual noises are observed. Do not use the instrument until it has been repaired by a Carl Zeiss Meditec authorised specialist.
- ❑ If, after switching on the instrument with the key switch, the display fails to light up, shut down the instrument, disconnect the power cable and attach a sign which clearly states that the instrument is out of order. Do not use the instrument until it has been repaired by a Carl Zeiss Meditec authorised specialist.
- ❑ Do not use force when inserting connectors. If the connector cannot readily be inserted into the jack or socket, make sure that you are using the correct connector for the socket. Do not pull the cables. Defective connectors should be repaired by a Carl Zeiss Meditec authorised specialist.

- ❑ Modifications and repairs to this instrument and any devices operated in combination with the VISULAS 532s are to be carried out by Carl Zeiss Meditec authorised specialists only. The manufacturer accepts no liability for damages caused by unauthorised access to the interior of the instrument. Such actions will also render any warranty claims invalid.
- ❑ This instrument may only be used with accessory equipment supplied by Carl Zeiss Meditec.
- ❑ For safety reasons and to ensure adherence to the guarantee conditions, the laser console may only be opened by specialist personnel with express written authorisation from Carl Zeiss Meditec.
- ❑ Do not use mobile phones or other devices not compliant with EMC-class B in the vicinity of this equipment, as their signals may cause the equipment to malfunction. The effect of radio signals on medical instruments is dependent on various factors and is therefore not predictable. To avoid electromagnetic interference, the instrument must be installed and operated as described in the operating instructions and using only those components supplied by Carl Zeiss Meditec.
- ❑ Disconnect the system from the power supply before opening the instrument, carrying out servicing or changing the fuses.
- ❑ Portable multiple sockets or extension cables should not be used.
- ❑ The manufacturer accepts no liability for damages arising as a result of unauthorised access to the interior of the instrument or the use of unauthorised accessories.
Such actions will also render any warranty claims invalid.

Warning

Internal components remain electrically live after the instrument has been switched off using the key switch.
To completely disconnect from the power supply, the instrument must be unplugged.

Requirements for operation of the instrument

The instrument will be installed by a Carl Zeiss Meditec authorised specialist. Please ensure that the following requirements for operation are always fulfilled.

- ❑ The specified environmental conditions for correct use are adhered to (see technical specifications)..
- ❑ The power cable is connected to a properly earthed power outlet.
- ❑ The power cable intended for use with this instrument is used.
- ❑ No external damage to the instrument is visible.

- ☐ All cables and connectors are in perfect working condition.
- ☐ Particular attention should be paid to warning and information labels, inscriptions and parts, such as screws or surfaces, which are marked in red.
- ☐ Do not cover ventilation grids or slits.

Always

- ☐ Use the lowest possible power necessary to achieve the desired effect.
- ☐ Focus the aiming beam as carefully as possible on the tissue to be treated.
- ☐ Never fire a laser pulse if you are unable to see the aiming beam clearly in the target area.
- ☐ Always switch off the aiming beam when the instrument is not in use.
- ☐ Because the aiming beam follows the same optical path through the beam delivery system as the therapy beam, it represents a good method for checking that the beam delivery system is undamaged. If the aiming beam spot at the distal end of the beam delivery system is not visible, is weak or appears diffuse, it may indicate damage to or improper function of the beam delivery system.

After using

- ☐ Always use the key switch to switch the instrument off.
- ☐ Always remove the key when the instrument is not in use.
- ☐ Do not allow unauthorised personnel access to the key switch key.

Safe working order

- ☐ This instrument is a high-quality technical product. To ensure optimum performance and that the instrument remains in safe working order, the instrument must be inspected by a Carl Zeiss Meditec authorised specialist at least once a year.

Warning

Optical fibres are not to be bent sharply, kinked or incorrectly fastened, as this may cause damage to the instrument and / or injury to the patient or user.

Laser safety

Safety precautions for the laser hazard area

This instrument is a class IV (4) laser device. The institution operating this instrument must take safety precautions to prevent any hazard. Applicable national and international regulations must be observed.

We draw your attention to some important aspects of these regulations:

- ❑ The laser hazard area is the area in which the maximum permissible exposure (MPE) can be exceeded. This includes areas into which the laser beam might unintentionally be deflected by a reflecting surface.
- ❑ The laser hazard area should be as small as possible and should be demarcated using suitable screens. Unauthorised persons should be prevented from entering the laser hazard area. The number of persons in the laser hazard area should be restricted to as few as possible. Where suitable screens are not in use, the entire treatment room must be considered as a laser hazard area.
- ❑ Warning signs must be placed at the entrance to the laser hazard area.
- ❑ A warning lamp must be installed at the entrance indicating when the laser is in operation.
- ❑ The laser hazard area must be clearly demarcated and marked as such during laser operation. A remote interlock contact can be fitted. This switches off the laser automatically if anyone enters the laser hazard area during operation. Further safety precautions should be clarified with local authorities or with your laser safety officer.
- ❑ All objects within the laser hazard area, including the floor, should have diffusely reflecting surfaces or should be covered with diffusely reflecting material.
- ❑ Only the patient and appropriately trained personnel may remain within the laser hazard area.
These persons must wear suitable laser safety goggles when the laser is in operation.
- ❑ Persons working in the laser hazard area must be given instruction concerning safety regulations, safety precautions and operation of the instrument at least once a year. A written record of this instruction, listing the participants, must be kept.

Safety officer

The institution operating the instrument and the safety officer are responsible for implementing all safety measures, such that the patient, the attending physician and any persons present during laser treatment are not exposed to any danger. In Germany, BGV B2 must be observed. In other countries, the relevant national regulations must be observed.

The responsibilities of the safety officer, who is to be appointed in writing by the institution operating the laser, include:

- ☐ Implementation of safety measures.
- ☐ Instruction of all persons working with the instrument in safety measures and correct operation of the instrument.
- ☐ Marking the laser hazard area.
- ☐ Checking the warning signs and equipment.
- ☐ Correct therapeutic use of the instrument.
- ☐ Safekeeping of the key switch key.
- ☐ Safekeeping of the laser device.
- ☐ Proper reconnection of the device after relocation.
- ☐ Correct keeping of the instrument logbook.

Laser safety goggles

- ☐ All persons present in the laser hazard area during laser operation must wear laser safety goggles.
The patient's eye which is not being treated must also be properly protected.
- ☐ Exceptions:
When observing the target area through
 - a surgical microscope equipped with a 532 nm physician's safety filter
 - a laser slit lamp supplied by Carl Zeiss Meditec
 - the VISULINK 532/U or
 - the LIO 532 laser indirect ophthalmoscope
 the physician does not require laser safety goggles.
- ☐ The institution operating the instrument is responsible for providing proper laser safety goggles. The safety goggles must be approved for the wavelength of the laser and for the type and intensity of the laser radiation (in Europe, in compliance with DIN EN 207).
Laser safety goggles for operation of the:
VISULAS 532s: D 450 to 532 L5

Note for VISULAS YAG III Combi

Laser safety goggles for operation in combination with the VISULAS YAG III: IRE 1064 L5

Patient safety

- ❑ The most important parameter in laser treatment is the power density at the target area, i.e. the applied laser power divided by the area of the laser spot. To achieve a constant physiological effect, the power output must also be changed when the spot diameter is changed. On the VISULAS 532s, the power output can be pre-set.

Note for VISULAS YAG III Combi

On the VISULAS YAG III Combi, the power output (VISULAS 532s mode) or the energy (VISULAS YAG III mode) can be pre-set.

Danger of explosion and fire

- ❑ Do not use the laser in combination with flammable anaesthetics.
- ❑ Keep explosive materials away from the laser hazard area. Easily flammable materials may cause fire.
- ❑ Laser beams can set fire to many explosive or flammable gases and liquids, including some solvents used in preparing for surgical procedures.
- ❑ Keep flammable drapes, surgical gowns, gauze and other flammable materials out of the beam path. We recommend using non-flammable materials and instruments and wearing gowns and clothing that do not easily catch fire.
- ❑ A fire extinguisher should be installed in the vicinity of the instrument.

Hazards caused by direct exposure and reflection

- ❑ This instrument emits visible laser radiation which is hazardous to the human eye and skin.

Note for VISULAS YAG III Combi

In combination mode, this instrument emits visible and invisible laser radiation which is hazardous to the human eye and skin.

- ❑ Do not expose the eye or skin to direct or scattered radiation! Metal and other surfaces reflect laser beams.
All persons present in the laser hazard area must therefore wear laser safety goggles as a precaution against unintentional direct or reflected laser radiation.
- ❑ Remove or cover any shiny, reflecting objects within the laser hazard area.
- ❑ Windows and reflecting walls should also be covered with non-flammable cloths.
- ❑ Safety precautions for dealing with noxious gases, dust and vapours, secondary radiation or explosive gas mixtures which may be produced by exposure of materials within the laser hazard area to laser radiation should be put in place.
- ❑ Only medical instruments which largely exclude hazardous reflections due to their design and surface finish should be brought into the laser beam path.

Faults

- ❑ If, during testing or operation of the laser system, faults occur which cannot be remedied by following the instructions given in the section *Maintenance, Miscellaneous*, the instrument should be switched off, the power cable disconnected and the instrument clearly labelled as being out of order.
Inform the responsible service technician.

Instrument safety systems

Safety equipment	Effect
Key switch (see 4, Fig. 14)	Unauthorised persons must be prevented from firing the laser. The user must remove the key after the treatment.
Remote interlock contact (see 9, Fig. 14)	The ability to fire the laser can be made dependent on the status of an external contact, for example on a door closer. A Carl Zeiss Meditec authorised specialist will be happy to assist, if you wish to connect the laser to a remote interlock. If the instrument is not connected to a remote interlock, a shorting plug which bridges the relevant contacts must be inserted in the interlock socket. The instrument is delivered with the shorting plug inserted.
Monitoring of therapy beam power	Laser firing is automatically locked when the metered power is outside the pre-defined range.
Safety monitoring	The electronic systems monitor various functions and values. If any of these values are outside the permitted range, the system prevents the laser from being fired.
LASER STOP button (see 3, Fig. 14)	The LASER STOP button is intended for control of unforeseeable hazardous situations. When this button is pressed, the system goes into Treatment / STANDBY mode and all functions on the control unit are locked. To resume operation, press the LASER STOP button again.
Laser warning lamp (see Fig. 29)	This lamp lights in Treatment / READY mode and in Treatment / STANDBY mode, when the aiming beam is on.

Warning and information labels

**Caution**

Observe all warning- and information labels!

Warning and information labels and their location on the instrument are shown in Fig. 1 to Fig. 4.

Should it become apparent that one of these labels is missing from your instrument, please contact Carl Zeiss Meditec or a Carl Zeiss Meditec authorised representative.


Fig. 1 Warning and information labels on the LSL 532s laser slit lamp


Fig. 2 Warning and information labels on the LSL YAG III Combi laser console


Fig. 3 Warning and information labels on the laser console


Fig. 4 Warning and information labels on the LIO 532 and VISULINK 532/U

Intended use

The VISULAS 532s can be used for the following applications:

- ☐ Photocoagulation of the retina
- ☐ Trabeculoplasty for treatment of glaucoma
- ☐ Iridotomy for treatment of glaucoma

Warning

Users in the USA:

The intended use is described on page 22 of the English language version of the operating instructions.

The VISULAS 532s system consists of the laser console as a radiation source, the LSL 532s laser slit lamp and the VISULINK 532/U and LIO 532 laser indirect ophthalmoscope as applicators. The laser console can also be used in combination with endo-probes.

The laser indirect ophthalmoscope can also be used for fundus examination. It is supplied together with a cold light source. Separate operating instructions from HEINE are included.

Warning

The VISULAS 532s may only be operated by physicians who have adequate knowledge of the medical application of the instrument, the effect of the laser beam on tissues and possible side effects.

Never fire the laser if you cannot clearly see the aiming beam spot in the target area.

Possible complications

Retinal bleeding and / or burns may occur, particularly where excessive laser power is used.

Note for VISULAS YAG III Combi:

The VISULAS YAG III Combi instrument system consists of the VISULAS 532s laser system, the VISULAS YAG III laser system and the LSL YAG III Combi laser slit lamp as a general applicator for the two laser systems.

For correct operation of the VISULAS YAG III, the VISULAS YAG III operating instructions must be observed.

Caution

For contraindications see the VISULAS YAG III operating instructions.

Contact lenses

Always use the correct contact lens for the intended application.

Warning

Always hold the contact lens at right angles to the laser beam. Incorrect positioning of the contact lens causes distortion of the laser focus. This can lead to unwanted side effects.

Never fire the laser if you cannot clearly see the aiming beam spot in the target area.

Operating principle

A laser crystal generates laser radiation, the frequency of which is doubled by means of a crystal.


Fig. 5 Generation of laser radiation with $\lambda = 532 \text{ nm}$

Radiation with a wavelength of 532 nm is strongly absorbed by bodily tissues. The thermal energy thus created causes coagulation of tissue protein, which can be used for medical purposes.

VISULAS 532s system design

The VISULAS 532s may consist of the following components:


Fig. 6 Components of the VISULAS 532s

Design of the VISULAS YAG III Combi

The VISULAS YAG III Combi may consist of the following components:


Fig. 7 Components of the VISULAS YAG III Combi

Controls, display elements and connectors

Laser slit lamp

- 1 Prism head, illumination and laser exit aperture**
- 2 Magnification changer**
- 3 Fixing screw**
This screw fastens the binocular tube or other component to the microscope body.
- 4 Interpupillary distance scale**
- 5 Binocular tube**
- 6 Eyepieces**
- 7 Adjustable eyecups**
The eyecups should be pushed in for spectacle wearers.
- 8 Laser spot size control knob**
- 9 Halogen lamp cover**
- 10 Slit length control knob**
1/3/5/9/14 mm
special slit positions (90°; $\pm 45^\circ$; 1 mm x 5 mm)
- 11 Slit width control knob**
(0 ... 14 mm)
- 12 Filter selector**
White line pointing up:
Full aperture combined with a heat-absorbing filter.
White line to the patient:
Blue filter for tonometry and fluorescence observations.
White line to the physician:
Green filter for examining the fundus.
- 13 Laser head YAG III**
- 14 Focus shift**
(--> *Focus shift and defocusing;*
see VISULAS YAG III operating instructions)
- 15 Hand switch for VISULAS YAG III**


Fig. 8 LSL 532s laser slit lamp


Fig. 9 LSL YAG III Combi laser slit lamp

Instrument base with headrest

1 Fixation lamp mounting

Caution

The laser slit lamp's optional fixation lamp is to be operated with the standard red blinking diode provided only. Otherwise, the fixation lamp may overheat and be damaged.

2 Patient headrest

3 Height adjustment for chinrest

4 Patient chinrest

5 Locking screw for instrument base

6 Laser slit lamp connection cable

Insert the cable into the socket (12, Fig. 14) on the laser console.

7 Mounting holes

The instrument base is mounted to the support using the appropriate screws.

8 Fast-action base lock

Lever tilted towards the patient: Base is fixed.

Lever tilted towards the physician: Base can be moved freely.

9 Brightness control for slit illumination

10 Joystick

X/Y fine adjustment by tilting the joystick;

(X/Y coarse adjustment is achieved by moving the complete instrument base).

11 Micromanipulator with locking ring

When the locking ring is unscrewed, the position of the slit image and laser spot can be finely adjusted. The VISULAS YAG III Combi has an additional hand switch in this position (see Fig. 9, 15). The VISULAS YAG III can be fired using a hand switch or foot switch.

12 Index mark for mid position

Indicates the height adjustment mid-position.

The height is adjusted by turning the joystick.

13 Eye level mark

14 Illumination swivel-angle scale


Fig. 10 Instrument base with headrest

Beam delivery in the laser slit lamp

The laser beam is interlaced into the illumination beam path, so that the slit image and the laser spot can move around the fundus in synchrony. This ensures that the area of the fundus which is to be coagulated is always illuminated, even in the peripheral fundus.

This means of combining illumination and laser beams is particularly advantageous for small illumination fields.


Fig. 11 Beam delivery in the LSL 532s laser slit lamp


Beam delivery in the YAG laser head (for VISULAS YAG III Combi)

A Nd:YAG laser (6) (mini-laser module) generates the therapy beam. The beam passes through the energy attenuator (2), is then combined with the aiming beam (1) and coupled into the viewing path of the corneal microscope by means of a coupling mirror (9).

The optical system (5) can be used to set the focus shift between the aiming beam and the therapy beam.

 Note:

The prism head of the slit illumination system must be lowered by about 30° before the YAG laser can be fired and the YAG laser beam can arrive at the patient's eye without encountering any obstacles.

- 1 4-point aiming beam diode laser
- 2 Energy attenuator
- 3 Energy meter and pulse counter
- 4 Partially transmitting mirror for energy meter
- 5 Focus shift lens system
- 6 Nd:YAG laser
- 7 Galilean magnification changer
- 8 Physician's safety filter
- 9 Mirror for coupling the laser beam into the viewing path
- 10 Slit lamp objective
- 11 Beam expander
- 12 Coupling mirror for the aiming beam

Fig. 12 Beam delivery in the YAG laser head

Laser console

The laser console (1) is the central power supply and control unit for the VISULAS 532s. Treatment parameters are set using the detachable control panel (2).

There is a fold-away handle (4) for carrying the console.

The key switch, the LASER STOP button, the optical fibre socket for the applicators, the PCMCIA card slot and the control panel connector are situated on the front panel of the console (see Fig. 14).

The applicator supply and control connectors and the power input are situated on the rear panel (see Fig. 14).



Fig. 13 Laser console with control panel

- 1 Laser console with control panel removed
- 2 Connector for control panel
- 3 LASER STOP button
- 4 Key switch
- 5 PCMCIA card slot
- 6 Optical fibre connector for applicators
- 7 Power supply for light source (not connected)
- 8 Foot switch connector
- 9 Connector for door interlock and external laser warning lamp
- 10 Serial port for servicing
- 11 Connector for YAG for combi laser use
- 12 Slit lamp connector
- 13 Power input connector

Fig. 14 Connectors and controls on the laser console


Fig. 15 Control panel, guide rod swivelled out

Control panel

The laser console is menu operated via the buttons on the control panel (1, Fig. 16) and the combined control knob / button (6).

The control panel can be attached to the laser console using the two guide rods (7 and 8) and the retaining elements (4 and 5) (snaps into place).

- To do so, first plug the connector (2) into the socket (3).
- Place the guide rods (7; 8) onto the retaining elements (4; 5) and push the control panel backward until it audibly clicks into place.
- To remove, hold the control panel at the front with both hands. Tilt the upper edge forwards so that the guide rods are released from the retaining elements.
- Remove the control panel by pulling it towards you carefully and disconnect the cable (2) by grasping the connector by the sleeve and pulling it directly backwards.

The lower guide rod (8) can be rotated by up to 90°, allowing the control panel to be safely angled towards the user (Fig. 15). This also allows the control panel to be inclined whilst mounted to the control panel.



Fig. 16 Removing the control panel

The combined control knob / button (6) is used to adjust the system parameters in the individual operating modes by turning the knob. Pressing the button selects the *laser power output* parameter in all modes.


VISULINK 532/U optical system

- 1 Optical fibre
- 2 Coupling lens
- 3 Beam expansion system
- 4 Objective
- 5 Reflecting mirror with physician's safety filter
- 6 Absorber
- 7 Fundus
- 8 Contact lens

Fig. 17 VISULINK 532/U optical system

VISULINK 532/U with accessories

- | | |
|--|--------------------------------|
| 1 Absorber | 10 Optical fibre |
| 2 Reflecting mirror with physician's safety filter | 11 SMA-connector |
| 3 Vertical mirror adjustment | 12 Cable clip |
| 4 Knurled ring for spot size adjustment | 13 Focusing aid |
| 5 Swivel arm | 14 Adapter for |
| 6 Spot size scale | Haag-Streit 900 BM / 900 BQ |
| 7 Spot size scale indicator mark | 15 Adapter for 20 / 30 SL |
| 8 Mount for mounting pivot | 16 Adapter for SL 115 Classic |
| 9 Swivel arm bearing pivot | 17 Adapter for SL 120 / SL 130 |

Fig. 18 VISULINK 532/U with accessories

LIO 532 optical system

1 Physician's eye**2 Lens +2 D**

The binocular body is fitted with +2 D lenses to assist the physician with accommodation. These lenses may be replaced with glass discs with a refractive power of 0 D.

3 Physician's safety filter**4 Binocular body**

The binocular body enables the physician to observe the fundus stereoscopically.

The laser beam arrives at the ophthalmoscope via an optical fibre and is directly coupled into the illumination beam path. The observation beam path (5), illumination beam path (7) and laser beam path (6) are then coincident.

5 Observation beam path**6 Laser beam path****7 Illumination beam path****8 Image plane and waist of the laser beam**

The physician observes the laterally-inverted image produced by the aspherical lens through the eyepieces of the binocular body.

The smallest laser spot diameter is achieved by selecting the working distance such that the focus ("waist") of the laser beam lies in the image plane.

9 Aspherical lens

The aspherical lens focuses the illumination light source onto the patient's iris.

If the 20 D aspherical lens is used, the total magnification of the system is approximately 2x.


10 Patient's eye

Fig. 19 Operating principle

LIO 532 laser indirect ophthalmoscope with accessories


Fig. 20 LIO 532 laser indirect ophthalmoscope and cold light source controls

- 1 Headband size adjuster**
People with smaller heads can replace the leather pad with the larger foam pad.
- 2 Clip for optical fibre** (on both sides)
The illumination optical fibre (**4**) is passed through one clip, the laser optical fibre (**6**) through the other.
The laser optical fibre contains anti-kink elements to protect it from excessive bending.
- 3 Headband height adjuster**
- 4 Illumination optical fibre** to HEINE cold light source
Please handle with care to avoid damage to the illumination optical fibre.
- 5 Sleeve nut** for the illumination optical fibre
- 6 Laser optical fibre** to VISULAS 532s laser console
Please handle with care to avoid damage to the laser optical fibre.
- 7 Eyepieces** with PD control
- 8 Height adjustment control for the illumination and laser beams**
- 9 Clamp for binocular body**
- 10 Mounting for binocular body**
- 11 Detent button** for binocular body mounting
- 12 Aspherical lens**
The laser indirect ophthalmoscope is supplied with a 20 D aspherical lens.
- 13 Power cable connector**
- 14 Cover**
The halogen lamp and a compartment for a spare lamp bulb are located under this cover.
Warning
Please observe the safety instructions in the operating instructions for the cold light source when opening the source or changing the halogen bulb.
- 15 Connection jack for illumination optical fibre**
The connector must click into place when the fibre is connected.
- 16 Brightness control**
- 17 Power switch**

Installation of the instrument

The instrument will be installed by a Carl Zeiss Meditec authorised specialist.

Caution

Optical fibres and application systems from other manufacturers are not to be used with this laser system.

- ☐ Avoid impacts and blows during transport of the instrument.
- ☐ Mount the system in a suitable location close to the physician and the slit lamp.

Caution

Ensure that the table or shelf on which the laser is placed is suitable for supporting a weight of 15 kg.

Caution

To avoid damage to the instrument, remove the key before carrying out this procedure.

- To install the laser on an IT 4L instrument table with lateral table mounting, the table mounting must first be screwed to the tabletop using the four screws and threaded holes provided.
- Now affix the laser slit lamp to the table.
- Slide the VISULAS YAG III (laser console), from the front, onto the upper shelf of the table mounting (for VISULAS YAG III Combi only).
- Finally, holding it vertically, slide the VISULAS 532s, from the front, onto the small side shelf.
- Place the foot switch in a suitable position under the table and insert the connector into the red socket (8, Fig. 14)


Note:

The cables and the corresponding sockets and ports are colour-coded for easy identification. The sockets and ports are also labelled.


Connectors and sockets are both marked with red dots. These must be lined up when connecting the cables.

- To disconnect the cable, grasp the connector by the sleeve and pull it directly backwards. Pulling at the sleeve releases the lock mechanism on the connector.

- 1 VISULAS 532s
- 2 Table mounting

Fig. 21 IT 4L with VISULAS 532s

- 1 Laser console
VISULAS YAG III:
- 2 VISULAS 532s laser console
- 3 Table mounting

Fig. 22 IT 4L with VISULAS YAG III Combi laser console

Caution

Only connect cables and sockets that are the same colour and size. Do not use force when inserting connectors.

- Install the remote interlock contact. If you do not wish to install an interlock, connect the yellow coding plug provided in place of the interlock (yellow connector, **9**, Fig. 14).
- The power cable for the IT 4L should be connected to a fixed power socket.
- Connect the console power cable to a fixed power socket or use one of the instrument table power sockets.

Note for VISULAS YAG III Combi:

- Connect the VISULAS YAG III to one of the sockets on the instrument table (see separate IT 4L operating instructions). The VISULAS 532s must be plugged into a separate fixed earthed socket.
- Remove the protective cap from the optical fibre of the applicator.

Caution

Never touch the end of the optical fibre with your fingers. Protect the end of the optical fibre from contamination.

Because of the laser's high power output, any dirt may lead to damage to the end surface of the fibre.

- Insert the optical fibre carefully into the socket for the optical fibre (**6**, Fig. 14) and screw in firmly.

Caution

Screw in by hand only. Do not use tools!

When operating with control panel removed:

- Fold out the control panel foot, unwind as much of the control panel cable as required and connect it to the front of the instrument (**2**, Fig. 14).
- Adjust the foot to the required angle and place the control panel in the location required.

When operating with control panel attached:

- Fold out the control panel foot, wind the control panel cable around the cable holder (**9**, Fig. 16) and attach the control panel to the laser console as described on p. 34.

Moving the IT 4L instrument table with the instrument

To avoid damage when moving the complete system, proceed as follows:

- Carefully lower the table to the lowest position.

Caution

Ensure that there are no objects in the path of the tabletop when moving the table downwards.

- The table can now be lifted slightly by the narrow side of the tabletop, on the slit lamp side, and placed onto the castors. Use both hands to lift the table.

Menu navigation

The VISULAS 532s is operated via a menu system with eight menu screens using the control panel.


Fig. 23 General menu navigation

VISULAS 532s description of menus

Note for VISULAS YAG III Combi

For a description of the menus for the VISULAS YAG III see the VISULAS YAG III operating instructions.

Power-up mode

Fig. 24 Power-up mode

- ❑ When the instrument is switched on using the key switch (4, Fig. 14) on the console, the start-up screen is displayed on the control panel (Fig. 24).
- ❑ The system automatically performs a self-test.
- ❑ The user interface includes a progress bar, which shows the progress of the system test.
- ❑ When the system test has been successfully completed, the program switches to the **main menu** (Fig. 26).
- ❑ You can switch directly into **Diagnosis mode** during the system test.

- ❑ To do so, touch the **eco** button


Diagnosis mode

This mode enables the user to use the LSL 532s laser slit lamp as a diagnostic slit lamp. The VISULAS 532s is used merely as a power source for the slit lamp.

- ❑ The laser module and temperature regulator remain switched off, meaning that the ventilators are operated at minimum speed (minimal noise production).
- ❑ The shutter remains closed in diagnosis mode.

Fig. 25 Diagnosis mode

- ❑ Touching the **ARROW-** button switches the system into **System test** mode.

Main menu

All **treatment mode** relevant parameters can be set in the main menu. These are:



Fig. 26 . Main menu

- ❑ Applicator selection (*slit lamp, LINK, LIO, ENDO probe*, for YAG laser in combination with LSL YAG III Combi).
- ❑ Program selection (three user-defined programs and a factory set default program for each applicator).
These programs can be used to set up pre-defined laser parameters (*power, pulse duration, repeat interval and aiming beam*).
- ❑ Activation / deactivation of the **audible signal** (VISULAS 532s only).
- ❑ Selection of **YAG firing mode** (hand / foot switch).
- ❑ The following modes can be accessed from the main menu.
 - Information mode
 - Diagnosis mode
 - Program definition mode
 - Treatment mode.
- ❑ Touch the **ARROW** button to switch to **Treatment / STANDBY** mode. The system will check whether the selected applicator is ready for operation.
If the applicator is not ready for operation, a message to this effect will be displayed.

The buttons in the **Main menu** have the following functions:

Function	Button	Function
Applicator selection		LSL 532s laser slit lamp
		VISULINK 532/U
		LIO 532
		Endo probe
		For YAG Laser in combination with LSL YAG III Combi
Program selection		Factory set default program
		P1, P2 and P3 are user-defined programs
Mode selection		Switch to information mode
		Switch to diagnosis mode
		Switch to program definition mode
		Switch to treatment mode
YAG firer		Laser is fired by hand (pushbutton on slit lamp joystick)
		Laser is fired using foot switch
Audible signal		Audible signal ON
		Audible signal OFF (VISULAS 532s only)

The reverse-coloured button (blue against white/yellow background) is activated.

System information

- ❑ Touch the **INFO** button to switch to **System information** mode.

This mode allows you to view system information. You can also set the *system time* (date and time), the *user interface language* and the *contrast* for the control panel LCD display.

The following system information is displayed:

- ❑ Company name / instrument name,
- ❑ Laser console and control panel software versions,
- ❑ Information on any faults which may have occurred,
- ❑ Operating hours counter / system time (time / date).

The following options can be set in **System Information** mode:

- ❑ Language. The user-interface is available in German, English, French and Spanish.
 - Touch the appropriate button to select the language required.
- ❑ LCD display contrast (to allow adjustment for the ambient brightness).
 - Touch the contrast button and adjust the LCD brightness display using the control knob (6, Fig. 16).

While **System information** mode is active, the shutter remains closed.

- ❑ Touching the **ARROW** button returns you to the **Main menu**.

Fig. 27 System information

english

Program definition mode

- ❑ Touch the **DEF** button to switch to **Program definition** mode.


Fig. 28 . Program definition mode

- ❑ In this mode, you can set the initial values for the laser parameters *power*, *pulse duration*, *interval* and *aiming beam intensity*. These initial values will be used whenever the system is restarted
- ❑ The software supports four programs.
Three user-definable programs (P1, P2, P3) are available for each of the applicators, *laser slit lamp*, *Link*, *LIO* and *ENDO*. They allow the user to define initial treatment parameter sets.
The program (P-) is a factory-set program with default values.
- ❑ To set individual laser parameters, touch the appropriate parameter button and select the value required by turning the control knob.
- ❑ Pressing the control knob takes you directly to the entry field for the laser power of the therapy beam. Rotate the knob to set the power required.

- ❑ Touching the **ARROW** button returns you to the **Main menu**.

Treatment mode

STANDBY

- ❑ Touching the **ARROW** button switches from the **Main menu** to **Treatment STANDBY** mode.

In treatment mode, the therapy beam can be applied via the selected applicator.

All treatment-related laser parameters are displayed on the control panel and can be changed in this mode.

- ❑ To set individual laser parameters, touch the appropriate parameter button and select the value required by turning the control knob.
- ❑ Pressing the control knob takes you directly to the entry field for the laser power of the therapy beam. Rotate the knob to set the power required.
- ❑ The pulse duration can be adjusted within the range 10 ms to 2500 ms or set to continuous ("cw"). The setting "cw" is only available for power outputs of 500 mW or less.

Warning

Pulse duration and interval are shown in ms (milliseconds). Older lasers show these values in sec (seconds).

100 ms is equivalent to 0.1 sec

The following restriction is placed on the interval setting in Repeat mode:

- An interval of between 100 ms and 6000 ms can be set, depending on power and pulse duration.


Fig. 29 . Treatment STANDBY mode; a) for LSL and VISULINK, b) for LIO and ENDO


- ❑ The aiming beam can be turned on and off by touching the **AIMING BEAM** button. When the aiming beam is on, this button is highlighted.

Aiming beam
ON

Aiming beam
OFF

- ❑ After touching the **BRIGHTNESS** button, the intensity of the aiming beam can be adjusted by rotating the control knob.
- ❑ A pulse counter is displayed. To reset the pulse counter to zero, touch the **RESET** button.
- ❑ Touching the **ARROW** button returns you to the **Main menu**.
- ❑ Touching the button for switching between treatment **STANDBY / READY** modes switches the system, after about three seconds, into **Treatment READY mode**. This transition time is visually indicated by a blinking button.

READY



Fig. 30 . Treatment READY mode; a) for LSL and VISULINK, b) for LIO and ENDO

- ❑ In **Treatment READY** mode, the aiming beam and laser warning lamp are always on.
- ❑ Parameters are set in the same way as for **Treatment STANDBY** mode.
- ❑ When the firing switch is activated, the laser beam is fired according to the parameters set.
- ❑ Touching the **STANDBY / READY** button returns you to **TREATMENT STANDBY** mode.
- ❑ Touching the **ARROW** button returns you to the **Main menu**.

Note:

If, in **Treatment READY** mode, more than 5 minutes elapses without the laser being fired, the system returns to **Treatment STANDBY** mode and emits a short audible signal.

Operation with laser slit lamp

Adjusting the slit lamp and eyepieces

Before treating a patient with laser radiation, you should be fully familiar with the use of the slit lamp. You should therefore carefully read the relevant user's manual, and in particular the safety notes.

- Ensure that the eyepieces are inserted fully and that the eyecups are pulled out (pushed in for spectacle wearers).
- Rotate the diopter setting ring on both eyepieces fully in the "+" direction (counter clockwise).
- Place a sheet of paper (business card or similar) in the object plane of the slit lamp.
- Centre the slit projector and binocular tube to the slit lamp base.
- Select the highest magnification setting on the slit lamp.
- Open the slit.
- Look through the eyepieces. Focus on the surface of the paper using the control lever and lock the slit lamp base.
- Close the illuminating slit to a narrow line.
- Select the lowest magnification setting on the slit lamp.
- Look through each eyepiece in turn, rotating the diopter setting ring for the eyepiece through which you are looking in the "-" direction (clockwise) until the surface of the paper is sharply in focus.
- Select the other magnification settings on the slit lamp.
- The image should remain sharply in focus at all magnifications. If it does not, the procedure should be repeated.
- Make a note of the eyepiece settings for future use. For future sessions, the eyepieces must simply be adjusted to this value.

**Note:**

If several physicians use this instrument, it is advisable to draw up a table with the individual refractive powers and keep it in a handy location near the instrument.

Positioning the contact lens

Hold the contact lens such that the laser beam meets the incident surface at a right angle. Incorrect positioning of the contact lens causes distortion of the laser focus. This can lead to unwanted side effects.

**Note:**

As a rule, you will get the best overall view if a magnification of 12x is selected on the slit lamp.

Laser treatment with VISULAS 532s and LSL

Positioning and focusing the laser spot

Exact positioning of the laser spot is crucial for successful laser treatment. This allows the laser power required to be minimised and significantly reduces unwanted side effects.

- Set the laser spot diameter by rotating the knurled ring (8, Fig. 8) on the slit lamp.
- Using the slit lamp joystick control, carefully position the aiming beam until the light spot appears sharply delineated and the tissue to be treated is seen sharply in focus.

Warning

Never fire the laser if the aiming beam is not unequivocally visible in the target area.

Note:

Always focus the slit lamp so that the retina and laser spot are both sharply in focus together, even if the laser spot appears smaller in other positions.

Laser configuration

- Connect the electric and optical connectors to the laser console.
- Switch on the VISULAS 532s laser console using the main power switch (4, Fig. 14).
- Set the laser parameters required in the menus. The menu screens and configuration procedures are described in the section "Description of menus" on p. 44.
- Operating the foot switch fires laser pulses using the parameters selected.

Note:

If no laser treatment is carried out, switch back to **Treatment / STANDBY** mode by touching the **STANDBY / READY** button. After each treatment session, switch off the instrument using the key switch and remove the key.

Operation with the VISULINK 532/U

The VISULAS 532s with VISULINK 532/U is compatible with the Carl Zeiss Meditec SL 120, SL 130, SL 115 Classic, 20 SL, SL 30 slit lamps and Haag-Streit slit lamp models 900[®] BM and 900[®] BQ.

Initial installation will be carried out by a Carl Zeiss Meditec authorised specialist.

Mounting the VISULINK 532/U on Carl Zeiss Meditec and Haag-Streit slit lamps

- The correct slit lamp adapter must be mounted on the slit lamp in order to attach the VISULINK 532/U.
- Insert the bearing pivot (9, Fig. 18) on the instrument's swivel arm (5, Fig. 18) into the hole in the slit lamp adapter (8, Fig. 18). To ensure easy mounting, ensure that the pivot is inserted into the guide mechanism as straight as possible. The VISULINK 532/U is engaged and locked into place, signalled by an audible click, through a slight pressure on the axle. The axle must be fully inserted into the bearing.
- The instrument can now be swivelled into the operating position. When swivelled into the operating position, the swivel mechanism will engage perceptibly, at which point it is locked in this position.

Warning


To avoid injury to the hands, do not grasp the swivel arm of the VISULINK 532/U too close to the axle when swivelling the unit in and out of position.

Note:

A swivelling absorber (1, Fig. 18) is fitted below the reflecting mirror on the VISULINK 532/U. This serves as a beam trap to remove unwanted radiation from the observation beam path. Please adjust the absorber so that visible reflections in the microscope are minimised and it is not possible for the beam trap to collide with the prism. Ideally the absorber should be horizontal.

- The VISULINK 532/U can be mounted at any angle of swivel.
- Connect the optical fibre to the laser console. Attach the cable clip to the vertical rod of the headrest. Ensure that the optical fibre is not kinked.
- The first time the instrument is used, the laser spot can, if required, be adjusted to the optical path of the slit lamp with the focusing aid provided. To do so, see pages 55 and 56 "Adjusting the slit lamp and VISULINK 532/U".

**Note:**

It is generally easier to mount and remove the VISULINK 532/U if the swivel arm is grasped close to the bearing pivot.

Removing the VISULINK 532/U from the slit lamp**Note:**

The VISULINK 532/U can be removed in the swivelled out position only.

- To do so, the instrument should be moved to the right until it is perceptibly at the limit stop. In this position, the swivel arm can be removed by lifting vertically out of the slit lamp adapter. An audible click should be heard.



**Caution**

When mounting and removing the VISULINK 532/U, ensure that the physician's safety filter does not collide with any part of the slit lamp.

**Note:**

When switching the VISULINK 532/U between different slit lamps, the slit lamp adapter should remain on the slit lamp. The link only should be changed. This makes repeated adjustment of the VISULINK 532/U unnecessary.

Fig. 31 Removing the VISULINK 532/U

- 1 Haag-Streit slit lamp
- 2 Slit lamp adapter
- 3 VISULINK 532/U
- 4 Slit lamp adapter
- 5 Zeiss SL 130

Fig. 32 VISULINK 532/U mounted on Haag-Streit slit lamp and Carl Zeiss Meditec SL 130

Adjusting the slit lamp and the VISULINK 532/U

The following procedure is used to ensure that the laser spot is focused in the object plane of the slit lamp.

Adjusting the slit lamp and eyepieces

- Ensure that the eyepieces are inserted fully and that the eyecups are pulled out (pushed in for spectacle wearers).
- Rotate the diopter setting ring on both eyepieces fully in the "+" direction (counter clockwise).
- Place a sheet of paper (business card or similar) in the object plane of the slit lamp.
- Centre the slit projector and binocular tube to the slit lamp base.
- Select the highest magnification setting on the slit lamp.
- Open the slit.
- Look through the eyepieces. Focus on the surface of the paper using the control lever and lock the slit lamp base.
- Close the illuminating slit to a narrow line.
- Select the lowest magnification setting on the slit lamp.
- Look through each eyepiece in turn, rotating the diopter setting ring for the eyepiece through which you are looking in the "-" direction (clockwise) until the surface of the paper is sharply in focus.
- Select the other magnification settings on the slit lamp.
- The image should remain sharply in focus at all magnifications. If it does not, the procedure should be repeated.
- Make a note of the eyepiece settings for future use. For future sessions, the eyepieces must simply be adjusted to this value.

Note:

If several physicians use this instrument, it is advisable to draw up a table with the individual refractive powers and keep it in a handy location near the instrument.

Adjusting the laser spot to the slit lamp

Note

After installation by a Carl Zeiss Meditec authorised specialist, the laser spot is already correctly focused for the image plane of the slit lamp. After switching to a different slit lamp, the laser spot can, if necessary, be readjusted to the slit lamp illumination using the focusing aid provided.


- First adjust the slit lamp and eyepieces as described above. Lock the base of the slit lamp.
- Select the highest magnification setting on the slit lamp.
- Activate the aiming beam and select a spot size of 100 μm .
- Attach the focusing aid to the mirror holder as illustrated in Fig. 34. Loosen the fixing screw (5, Fig. 33) of the adapter by about half a turn, then move the adapter and VISULINK 532/U until the white structured lines are visible in sharp focus. Tighten the fixing screw. Recheck the focus and repeat the adjustment process if necessary.
- The VISULINK 532/U is now adjusted to this slit lamp and does **not** require readjustment as long as the adapter is not removed from this slit lamp.
- Centre the laser spot to the slit:

Vertically:

- Adjust the illumination slit so that it is horizontal.
- The laser spot is adjusted to the centre of the illumination slit using the dial (1, Fig. 33).

Horizontally:

- Adjust the illumination slit so that it is vertical.
- The laser spot is adjusted to the centre of the illumination slit by rotating the mirror holder (3, Fig. 33) in its entirety.


- 1 Dial for vertical adjustment
- 2 Reflecting mirror
- 3 Mirror holder
- 4 Adjusting ring with scale for spot size adjustment
- 5 Fixing screw

Fig. 33 Laser spot adjustment

Warning

Adjustment using equipment other than that or procedures other than those described here may result in exposure to hazardous radiation.

- 1 Focusing aid
- 2 Projection surface
- 3 Bearing bracket
- 4 Stop pin

Fig. 34 Mounting the focusing aid on the VISULINK 532/U

Laser treatment using the VISULINK 532/U

Before carrying out laser treatment, you should be thoroughly familiar with the operation of the slit lamp as a diagnostic instrument and with the VISULINK 532/U and VISULAS 532s.

- Start the laser console as described on p. 40.
- Insert the laser optical fibre into the light optical fibre socket (6, Fig. 14) on the laser console and screw in firmly.

Caution

Do not use tools to screw in the optical fibre!

Warning

Check the integrity of the reflecting mirror with physician's safety filter (2, Fig. 18) before each use.

Significant contamination or detached coatings can lead to exposure to hazardous radiation.

- Select the LINK applicator from the control panel.
- Select the appropriate program with the required operating parameters, as described on p. 48, or set the parameters in **Treatment / STANDBY** mode.
- Touch the **ARROW** button to switch to **Treatment / STANDBY** mode.
- Confirm the safety question by touching the **TICK**.

Caution

Never begin treatment if a different applicator is displayed.

- Select the laser spot diameter required by rotating the knurled ring (6, Fig. 18) on the VISULINK 532/U.
- Switch the aiming beam on.
- Hold the contact lens such that the laser beam meets the incident surface vertically. Incorrect positioning of the contact lens causes distortion of the laser focus. This can lead to unwanted side effects.

Note:

Because of the small working distance, it may not be possible to use the longer versions of some contact lens. Please use the shorter versions of these (e.g. Goldmann three mirror-contact lens, short, see p. 81).

- Position the aiming beam until the light spot appears sharply delineated and the tissue to be treated is seen sharply in focus.
- Touch the **STANDBY / READY** button to switch to **Treatment / READY** mode.
- Recheck the treatment site and operating parameters.
- Fire the laser pulse by operating the foot switch.

**Note**

In rare cases, you may perceive chain-like multiple reflections of the aiming beam. These are virtual reflections, i.e. they are artifacts of the microscope only and are not present on the retina. The actual aiming beam is easily distinguishable by its brightness. These multiple treatment laser images do not represent a risk to the eyes of the physician carrying out the treatment or the patient.

**Warning**

Never activate the laser if irregular multiple reflections of the aiming beam are visible.

Never fire the laser if the aiming beam is not unequivocally visible in the target area.

- After treatment, switch off the laser console using the key switch and remove the key.

Operation with the LIO 532 laser indirect ophthalmoscope

Before carrying out laser treatment, you should be thoroughly familiar with the operation of the laser indirect ophthalmoscope as a diagnostic instrument.

To prevent dust from entering the optical system, do not remove the laser optical fibre from the laser indirect ophthalmoscope.

For safety reasons, do not switch on the laser console unless you intend to perform laser treatment.

Connecting the cold light source

- Set up the cold light source in a suitable position to ensure unobstructed air circulation.
- Insert the illumination optical fibre (4, Fig. 20) into the connector (15, Fig. 20) on the cold light source. The connector must click into place.
- Switch on the cold light source and set the brightness to maximum.

**Warning**

Do not look directly into the light source (connector on the source or illumination optical fibre)!

The cold light source must be plugged into a separate fixed earthed socket. (Portable multiple sockets must not be used).

Adjusting the head strap

- Adjust the eyepieces to your approximate interpupillary distance. The eyepieces should be symmetrically positioned about the axis of the PD scale.

- Put on the laser indirect ophthalmoscope and adjust the head strap using the size (1, Fig. 20) and height (3, Fig. 20) adjusters. The device should sit firmly on the physician's head. People with smaller heads can replace the leather pad with the larger foam pad.
- Loosen the clamp (9, Fig. 20) and adjust the height and tilt of the binocular body, such that the eyepieces are as close as possible to the eyes and allow comfortable viewing.
Observe the light spot on a surface at a distance of about 30 to 40 cm.
- Centre the light spot vertically using the control knob (8, Fig. 20). Centre the light spot horizontally by correct the positioning of the headset if required. Operators who wear bifocal spectacles should look through the far vision segment of the spectacles.

Adjusting the binocular body

- Adjust the interpupillary distance by observing the light spot alternately through the left and right eyes and adjusting the eyepieces (7, Fig. 20) until the light spot appears precisely at the centre of the field of view.
- Take off the laser indirect ophthalmoscope and adjust the eyepieces so that the interpupillary distance remains the same and the eyepieces are arranged symmetrically around the PD scale. Make a note of this value to facilitate future adjustment.
- Put the laser indirect ophthalmoscope back on and check that the PD setting has not changed.
- After the head strap and binocular body have been properly adjusted, depress the detent button (11, Fig. 20) and swing the binocular body upward. When you swing the binocular body back down, it should return to the same working position. If it does not, the head strap is too loose. Correct the positioning using the adjusters (1, Fig. 20) and (3, Fig. 20).

Examining the fundus

- Prepare the patient, for instance, by dilating the pupils.
- Set the illumination to a low intensity. Always start with a low light intensity to protect the patient's eyes and increase lamp life.
- Hold the aspherical lens (12, Fig. 20) in front of the patient's eye so that the engraving on the lens mount faces the operator. Let your little finger rest on the patient's head to support your hand. The distance between operator and patient should be about 40 cm.
- Centre the light spot in the patient's pupil. Bring the image of the fundus into focus by moving your head forward or backward and changing the position of the aspherical lens.


Fig. 35 Replacing the lenses

Replacing the +2 D lenses

The binocular body is fitted with +2 D lenses to assist the physician with accommodation. If required, these lenses may be replaced with the 0 D glass plates provided.

- Remove the eyecups.
- Unscrew the eyepieces.
- Lift out the +2 D lenses by the lateral slots using a pointed instrument. Ensure that you do not scratch the lenses and that the physician's safety filter, which is situated directly below the lens, is not loosened or damaged.
- Insert the 0 D glass plates.
- Screw the eyepieces back on and reattach the eyecups.

Laser treatment with the LIO 532

Before carrying out laser treatment, you should be thoroughly familiar with the operation of the laser indirect ophthalmoscope as a diagnostic instrument and with the VISULAS 532s laser console.

Preparation

- Start the laser console as described on p. 40.
- Insert the laser optical fibre into the light optical fibre socket (6, Fig. 14) on the laser console and screw in firmly.

Caution

Do not use tools to screw in the laser optical fibre!

Warning

Check the integrity of the lenses and glass surfaces before each use. Significant contamination or detached coatings can lead to exposure to hazardous radiation.

- Prepare the patient for the treatment.
- Turn on the laser console using the key switch.
- Connect the laser indirect ophthalmoscope to the cold light source and turn on the source.
- Put on and adjust the laser indirect ophthalmoscope.

Laser treatment

- Select the **LIO 532** applicator from the control panel.
- Select the appropriate program with the required operating parameters, as described on p. 48, or set the parameters in **Treatment / STANDBY** mode.
- Activate or deactivate the audible signal.
- Touch the **ARROW** button to switch to **Treatment / STANDBY** mode.
- Confirm the safety question by touching the **TICK**.

Warning

Never begin treatment if a different applicator is displayed.

- Switch on the aiming beam.
- Observe the fundus of the patient's eye through the indirect ophthalmoscope. The aiming beam should be visible as a red spot in the centre of the field of view.
- Move your head forwards and backwards until you obtain the smallest spot size.

The desired physiological effect may fail to be achieved if the laser spot is not properly focused.

If the 20 D aspherical lens is used, the diameter of the laser spot will be 360 µm (for an emmetropic patient eye).

- Touch the **STANDBY / READY** button to switch to **Treatment / READY** mode.
- Recheck the treatment site and operating parameters.
- Operate the foot switch to fire the laser pulse.

Warning

Never fire the laser if the aiming beam is not clearly visible in the target area.

- After treatment, switch off the laser console using the key switch and remove the key.

Warning

If, after using the indirect ophthalmoscope, you change to a different applicator (e. g. slit lamp), the correct system must be selected from the main menu.

Operation with the intraocular probe (Endoprobe)

Description

Each intraocular probe is supplied in a sterile pack. The probe consists of the components illustrated in Fig. 36.


Fig. 36 Intraocular probe

Warning

The intraocular probe is supplied in a sterile condition. To preserve sterility, the following measures must be observed:

- ☐ Observe storage requirements and expiry dates.
- ☐ Carefully store the sterile probe pack in the cardboard box provided. Remove it immediately before use.
- ☐ Do not use if the pack or probe is damaged!
- ☐ The sterile pack should be opened immediately before use.
- ☐ The probe is for single use only and should be disposed of immediately after use. Resterilisation is not permitted. Reuse of probes may cause infection.

When disposing of probes, the regulations for handling infectious material must be observed.

Optical fibre fatigue as a result of resterilisation and reuse may result in fibre breaks and misdirection of the laser beam. This may result in injury.

- ☐ Connection of the intraocular probe to the VISULAS 532s laser console is described in section *Application* on p.64.

Physician's safety filter

Warning

To protect the attending physician's eyes from laser radiation, the surgical microscope must be equipped with a physician's safety filter.

Fixed physician's safety filters (000000-1272-301), which remain in the beam path throughout the treatment and active physician's safety filters (000000-1336-953, simple and 000000-1338-594, double), which move into place when the laser is fired only, are available.

The physician's safety filter should be installed by a trained and authorised service technician.

Please observe the following points:

- ❑ The physician's safety filter (for 532 nm) protects the attending physician's eyes when looking through the binocular tube of the surgical microscope.
- ❑ If an assistant's scope is connected, the physician's safety filter must be mounted upstream of the beam splitter to the assistant's scope.
- ❑ On OPMI CS surgical microscopes, the beam splitter is built into the body of the microscope. It is therefore not possible to mount the physician's safety filter upstream of the beam splitter. If you wish to mount an assistant's scope on the beam splitter, you must therefore install a second physician's safety filter between the assistant's scope and the straight binocular tube.
- ❑ If the surgical microscope is equipped with an assistant's microscope, a second physician's safety filter must be built into the optical path of the assistant's microscope.


Fig. 37 Intraocular probe

Application

- Open the pack and carefully remove the SMA connector of the sterile intraocular probe. The handpiece can be left in the pack until it is needed.
- Remove the protective cap (2, Fig. 36) from the optical fibre adapter.
- Screw the SMA connector (1, Fig. 36) into the laser console's optical fibre socket (6, Fig. 14).

Caution: Screw in the SMA connector by hand only. Do not use tools!

- If you need to change the intraocular probe during treatment, the laser console must be switched to **Treatment / STANDBY** mode.

Sterile operation

To set the laser parameters under sterile conditions, the control panel can be operated using a sterile swab or a sterile stylus (3, Fig. 38).

Caution


Do not use sharp objects or pens to operate the control panel.

To maintain the sterility of the sterile end of the swab or sterile stylus between uses, a holder (2, Fig. 38) can be attached to the control panel. Insert the sterile stylus into the holder when it is not required.

Warning

The swab or sterile stylus should be held by the sterile end only. Ensure that the sterile end does not come into contact with non-sterile materials or persons. Dispose of the swab or sterile stylus properly after use.

- 1 Control panel
- 2 Sterile stylus holder
- 3 Sterile stylus

Fig. 38 Control panel with sterile stylus and holder

Warning

The sterile stylus is supplied in a sterile state together with the intraocular probe. To preserve sterility, the following measures must be observed:

- ☐ Observe storage requirements and expiry dates.
 - ☐ Carefully store the sterile probe pack in the cardboard box provided. Remove it immediately before use.
 - ☐ Do not use if the pack is damaged!
 - ☐ The sterile pack should be opened immediately before use.
 - ☐ The sterile stylus is for single use only and should be disposed of immediately after use. Resterilisation is not permitted. Reuse of sterile styli may cause infection.
- When disposing of sterile styli, the regulations for handling infectious material must be observed.

Laser treatment

Warning

Before treatment check that the physician's safety filter is fitted in the surgical microscope.

The physician's safety filter is labelled with a continuous green strip and the words "VISULAS Safety Filter 532 nm".

- Start the laser console as described on p. 40.
- Switch on the VISULAS 532s using the key switch.
- Select the ENDO-probe as applicator and touch the **ARROW** button.
- Select the physician's safety filter used "**Fixed**" (1, Fig. 37) or "**Active**" (2 or 3, Fig. 37).


Fig. 39 Screen for confirming that physician's safety filter is fitted

- Select the appropriate program with the required operating parameters, as described on p. 48, or set the parameters in **Treatment / STANDBY** mode.

Warning

Never begin treatment if a different applicator is displayed or if there is no physician's safety filter fitted.

- Connect the intraocular probe to the laser console as described in *Application on p. 64*.
- Switch on the aiming beam.
- Hold the tip of the intraocular probe near a surface and adjust the brightness of the aiming beam.
- Introduce the intraocular probe.

- Switch to **Treatment / READY** mode. The foot switch is now activated.

Warning

Only operate the foot switch, if you have clearly localised the target area and are sure that you wish to use the therapeutic laser beam.

- Localise the target area and adjust the size of the laser spot.
- The tip of the probe should normally be 2 to 5 mm from the retina, depending on the spot size and power density required.

Warning

Do not fire the laser beam if the tip of the intraocular probe is touching the tissue.

- Operating the foot switch fires treatment laser pulses using the parameters selected.
- After treatment, switch off the laser console using the key switch and remove the key.

Troubleshooting tables

Faults are indicated as system information on the control unit display. The user is also alerted by short audible signals.

Note: Follow the instructions given.

Warning message: Remedy the fault according to the troubleshooting tables below.

If a fault occurs, which cannot be remedied by using the following tables, please contact a Carl Zeiss authorised specialist, citing the error number displayed.

Label the instrument as being out of order and remove the key from the key switch.

Warning messages and instructions

Message	Remedy
The interlock is activated, please close the door.	Shut the door or check that the shorting plug is correctly plugged into the yellow socket.
The LASER STOP is activated; to continue, please unlock.	Unlock the LASER STOP button.
Slit lamp or physician’s safety filter not connected!	Check that the slit lamp and the surgical microscope (ENDO physician’s safety filter) are correctly connected and properly inserted.
Optical fibre not connected!	Check that the optical cable is properly inserted and screwed in tightly.
Foot switch not connected.	Check that the foot switch is correctly connected and that the cable is properly inserted.
Temperature exceeded. Please wait!	If this fault occurs frequently, ensure that the instrument’s ventilation slots are not covered and that the ambient temperature falls below 35 °C.
Adjusting pulse interval.	If the laser parameters ‘pulse duration’ or ‘power’ are changed, the system may carry out an adjustment of the ‘pulse interval’ parameter. To avoid overloading the system, the interval between laser pulses is increased.
Power deviation exceeds 20 %!	The internal monitoring systems have identified that the actual laser power output deviates from the target power output selected on the control panel by more than 20%. Should this error message occur repeatedly, please inform a Carl Zeiss Meditec authorised specialist.
Release the firing switch.	Release the foot switch.
Monitoring systems deactivated. For service use only!	You are in a special service mode. All error messages are deactivated. Please contact our service department. Patients must not be treated in this mode!
LINK applicator selected. Do you want to use LINK?	Safety question, to check that the correct applicator has been selected. Confirm by touching the <i>tick</i> , or abort by touching the <i>‘X’</i> .

Message	Remedy
LIO applicator selected. Do you want to use LIO?	Safety question, to check that the correct applicator has been selected. Confirm by touching the <i>tick</i> , or abort by touching the 'X'.
Please select the physician's safety filter used and check that it is fitted.	Safety question, to check that the physician's safety filter is fitted in the surgical microscope. If you are using a fixed safety filter, touch 'fixed safety filter 532 nm'. If you are using an active safety filter, touch 'active safety filter 532 nm'. If no safety filter is fitted, abort by touching the 'X'.
Error 15	A button has been pressed for more than 5 seconds. Release the button and wait for a moment.
In the STANDBY and READY modes, demo is displayed beside the aiming beam button. The laser cannot be fired.	Switch to STANDBY mode. Touch the demo button until it disappears (about 5 seconds).

Note for VISULAS YAG III Combi

Swivel lamp casing out of position	Swivel the lamp casing out of position. Check that all the connectors on the rear of both consoles are inserted properly.
In the STANDBY and READY modes, during YAG operation, demo is displayed beside the aiming beam button. The laser cannot be fired	Switch the VISULAS YAG III off and then on again. While the progress bar is active during the YAG III system test, touch the ZEISS logo on the screen, then press the control knob/button.

Faults during operation of the LIO 532

Fault	Possible cause	Remedy
No or weak illumination	- Cold light source power cable not connected. - Cold light source not switched on. - Brightness too low. - Power failure - Halogen lamp in cold light source defective. - Illuminating optical fibre defective.	- Connect cold light source power cable. - Switch on cold light source. - Adjust brightness. - Call in-house electrician. - Replace halogen bulb. - Replace illumination optical fibre.
Strong reflections visible in the image.	- Aspherical lens held the wrong way.	- Hold aspherical lens so that its engraving faces the operator.
No laser beam.	- Laser console in Treatment STANDBY mode. - Error message or demo mode on laser console. - Laser optical fibre defective. - Laser coupler defective. - LASER STOP pressed.	- Switch to Treatment READY mode. - Refer to VISULAS 532s laser console trouble-shooting table. - Contact the service department. - Contact the service department. - Release the LASER STOP button.

Faults during operation of the laser slit lamp

Fault	Possible cause	Remedy
No function at all	- Laser console power cable not connected. - Connecting cable not properly inserted into laser console. - Fuse defective. - Control unit not switched on.	- Connect power cable. - Insert connector into socket (12, Fig. 14). - Replace fuse (see also VISULAS YAG III operating instructions). - Switch on the control unit.
No slit illumination	- Slit closed. - Halogen bulb defective.	- Adjust slit width using knob (11, Fig. 8). - Replace halogen bulb (--> <i>page 73</i>).
Poor slit illumination	- Selected lamp voltage too low. - Halogen bulb not properly inserted. - Slit length not correctly adjusted. - Slit width not correctly adjusted. - Filter selector not correctly set.	- Set halogen lamp brightness as required using the brightness control (9, Fig. 10). - Ensure that halogen bulb is properly inserted (--> <i>page 73</i>). - Adjust slit length using knob (10, Fig. 8). - Adjust slit width using knob (11, Fig. 8). - Check position of filter selector knob (12, Fig. 8) (white line should be pointing upward).
Difficulties when looking through corneal microscope	- Corneal microscope magnification not correctly adjusted. - Interpupillary distance of binocular tube not correctly adjusted. - Eyepieces not correctly adjusted.	- Turn magnification changer knob (2, Fig. 8) to select the desired magnification. - Adjust interpupillary distance of binocular tube. - Adjust eyepieces (--> <i>page 51</i>).
Movement of instrument base too stiff	- Instrument base locked with fast-action base lock. - Base locked with locking screw.	- Unlock fast-action base lock (8, Fig. 10). - Unscrew locking screw (5, Fig. 10).

**Note for VISULAS YAG III Combi**

If, in YAG mode, an error occurs on the VISULAS YAG III, an error message will be displayed on the control panel.
At the same time, an error message with a specific error code will be displayed on the YAG console.
(see also VISULAS YAG III operating instructions).

An error has occurred on the VISULAS YAG

Faults during operation of the VISULINK 532/U

Fault	Possible cause	Remedy
No aiming beam visible.	- Aiming beam not activated. - Link swivelled out. - Link not connected to laser console. - The intensity of the slit lamp is too high.	- Switch on the aiming beam. - Swivel the link into position. - Connect the end of the optical fibre to the laser console's output socket. - Reduce the intensity of the slit lamp to the required level.
Aiming beam is not at the centre of the illumination slit.	- Reflecting mirror is not correctly adjusted to the illumination slit.	- Adjust the reflection plate to the illumination slit as described on p. 56.
A precisely repeating series of aiming beams is visible.	- Virtual reflections, visible with the microscope only, have been produced	- The therapeutic beam is not affected. Read the instructions at the top of p. 58 and continue working accordingly
Irregular multiple reflections of the aiming beam are visible	- The reflecting mirror on the Link is contaminated. - The coating on the reflecting mirror is detaching.	- Clean the reflecting mirror as described on p. 74, <i>Cleaning of optical components</i>. - Stop using the equipment immediately and call a Carl Zeiss Meditec authorised specialist. Mark the equipment as "out of order". Do not carry out any further treatments.
The field of view is partly overlaid with scattered light.	- The absorber plate at the lower end of the link is contaminated. - The angle of inclination of the absorber plate is not correctly adjusted for the slit lamp.	- Clean the absorber plate as described on p. 74, <i>Cleaning of optical components</i>. - Swivel the absorber plate into position so that no scattered light effects are visible.
You are dazzled by bright treatment light in the slit lamp microscope during treatment.	- Physician's safety filter in reflecting mirror faulty.	- Stop treatment immediately and contact a Carl Zeiss Meditec authorised specialist as soon as possible.
Difficulty in using the contact lens.	- Contact lens is too long.	- Use a shorter context lens.
Link cannot be used.	- There is a foreign body in the quick swap station. - Link is misaligned.	- Remove the foreign body. - Insert the pivot vertically into the quick swap station.
Link cannot be removed.	- Link is not in the correct position. - Link becomes misaligned when trying to remove it.	- Swivel the Link out completely until it clicks perceptibly into place. - Grasp the Link close to the pivot and remove it vertically from the mounting.


Changing the halogen bulb in the laser slit lamp

Warning

Before changing the bulb, switch off the instrument using the laser console key switch and disconnect the power cable. Allow the lamp to cool down or wear protective gloves.

- Replace the halogen bulb (3) as illustrated in Fig. 40.
- To unlock the lamp holder, swing lever (5) aside.
- Do not touch the bulb with bare fingers.
- When inserting the bulb, make sure the lamp holder (2) is properly positioned in relation to the locating pins (4).

- 1 Halogen lamp cover
- 2 Lamp holder
- 3 Halogen bulb 12 V / 30 W
- 4 Locating pins
- 5 Lever

Fig. 40 Changing the halogen bulb

Care of the instrument

Warning

Disconnect the power cable before cleaning.

Cleaning of optical components

The multilayer T* (T star) antireflection coating of the optical components (e.g. eyepieces and objectives) results in excellent image quality. Even slight contamination of the optics or a fingerprint will impair image quality. In order to protect internal optical elements from dust, never store the instrument with objective, binocular tube or eyepieces removed. After use, the instrument should be covered to protect it from dust. Always store objectives, eyepieces and accessories that are not being used in dust-free cases.

The external surfaces of optical components (eyepieces, objectives), can be cleaned as required:

- Dust on the optical surfaces can be removed using a squeeze blower or a clean, grease free brush.
- Thorough cleaning can be performed quickly and easily using a moist, anti-static cleaning cloth.
Please follow the instructions on the cleaning cloth packaging.

This will enable you to maintain the high image quality of your product.

Cleaning and disinfection of painted surfaces

All painted surfaces can be cleaned with a damp cloth.

Do not use aggressive or abrasive cleaning agents.

Use spray or wipe disinfectants for cleaning and disinfecting the instrument casing, foot switch, chinrest and forehead strap.

The display should be cleaned with a lightly moistened cloth only.

Note:

Ethanol or isopropanol can be used to clean the control panel.

Ensure that no moisture penetrates the instrument during cleaning and disinfection.

Paper pads should be replaced after each patient.

Replacing the LIO 532s illumination optical fibre

- Switch off the cold light source and disconnect the optical fibre from the source.
- To disconnect the optical fibre (4, Fig. 20) from the binocular body, unscrew the sleeve nut (5, Fig. 20).
- Connect the replacement optical fibre in the same way. Ensure that the thread is not cross-threaded when screwing in the sleeve nut.

Cleaning the pads of the LIO 532s headstrap

The forehead and back of head pads can be removed and cleaned with a mild soap solution.

Storing the indirect ophthalmoscope

When not in use, to prevent optical elements from being contaminated by dust, always store the laser indirect ophthalmoscope in its carrying case.

- Insert the illumination optical fibre connector into the holes provided in the foam rubber.
- Screw the protective cap onto the end of the laser optical fibre.
- Adjust the headband height and size to minimum.
- Undo the clamp (9, Fig. 20) for the binocular body and place the laser indirect ophthalmoscope in the case. The binocular body should be placed as far to the rear as possible to prevent the laser optical fibre from being damaged when the case is shut.

Caution

When shutting the case, ensure that the optical fibre does not become caught or kinked.

Transportation instructions

Before transporting the instrument, disconnect the individual components by disconnecting the connecting cables.

Use an appropriate means of transport, such as the transport case, in order to prevent hard impacts and mechanical damage.


Reconnect the cables at the new installation site

(--> *Controls, display elements, connectors, p. 26 ff., p. 40*).

Disposal

This instrument contains electronic components. At the end of its lifetime, the instrument should be disposed of in accordance with the relevant national regulations.

Disposal of the product within the EU

EU directives and national regulations currently in force at the time of marketing prohibit the disposal of the product specified on the delivery note in domestic waste or by municipal waste disposal companies.

For more information on disposal of the product, please contact your local dealer, the manufacturer or their legal successors. Please also take note of the manufacturer's current online information.

If the product or its components are resold, the seller has the duty to notify the buyer that the product must be disposed of in accordance with currently valid national regulations.

Safety inspections

To ensure proper functioning of the instrument, it should be subjected to regular safety inspections.

These safety inspections should be performed annually by the manufacturer or a manufacturer authorised specialist and the results recorded in the instrument logbook.

These inspections should also be carried out when the instrument has been out of use for more than one year.

The following points should be subject to inspection / the following procedures carried out:

Measure	Procedure
Electrical safety	Protective earth resistance with power cable: $R_{PE} < 0.2 \text{ ohm}$. Measurements should be carried out on all conductive parts which could carry hazardous voltages in the event of a fault.
	Earth leakage current: Earth leakage current is to be measured while the device is in operation using a standard meter (Measurement instructions must be observed!). 0.5 mA in acc. with IEC 60601-1: $\leq$ or (in Germany only) 1 mA stand-by unit leakage current in acc. with DIN VDE 0751: $\leq 1 \text{ mA}$.
	Visual inspection in accordance with the checklist "Inspection of medical devices" (safety folder).
Safety inspection	Foot switch
	Visual inspection of physician's safety filter
	Safety switch-off, if actual power deviates from nominal power.
	Remote interlock contact
	Function test
Acceptance protocol	Check adjustment of the therapy laser to the slit lamp.
	Check aiming beam and therapy laser beam power downstream of the applicator.
	Check the internal power meter.
Cleaning the slit lamp lens	
Preparation of the instrument for use and handing it over to the user.	

Note for VISULAS YAG III Combi:

Safety inspection	Safety switch-off, if the actual energy and power deviate from the nominal energy and power
	Threshold values for YAG pulse modes 1, 2, 3
	Optical breakdown in air
	Hand switch

Calibrating the power measurement system

The VISULAS 532s is calibrated such that the power setting on the control panel coincides with the actual power delivered to the target area. Calibration of the power measurement system must be checked at least once every 12 months. This is performed as a routine measure by our service technicians during servicing. The US Food and Drug Administration (FDA) requires that manufacturers of Class III and IV medical laser instruments provide US customers with power measurement system calibration procedures.

**Warning**

Power measurement system calibration is a service procedure to be performed by trained and authorised service technicians only. Because of the risk of exposure to laser radiation, this procedure is dangerous. Improper calibration can lead to serious injury.

The control unit and laser console may only be opened by Carl Zeiss Meditec service technicians or by other persons who have been expressly authorised to do so in writing.

The manufacturer accepts no liability for damages or injury caused by unauthorised access to the interior of the instrument. Such actions will also render any warranty claims invalid.

These instructions are written for service technicians who have been officially trained by the manufacturer. Possession of service literature (service manual, installation instructions or similar) and/or the required tools does not entitle unauthorised persons to calibrate or otherwise service this laser device.

When properly calibrated, the system will compensate for any power losses in the application system used. It is therefore critical that the application system is adjusted in accordance with the regulations.

Calibration procedure

Tools and aids required for calibration of the laser system:

- ☐ Coherent LM10 calibrated power measurement system for measurement of radiation at 532 nm within a range of 0.05 - 10 W
- ☐ Dummy plug with adapter
- ☐ Laser safety goggles for 532 nm laser radiation
- ☐ Service plug for serial port

Proceed as follows:

Warning

All persons in the laser hazard area while the laser is in operation must wear laser safety goggles.

Calibrating the laser console

- Disconnect the optical fibre from the laser console.
- Connect the measuring head using the adapter.
- Connect the adapter to the optical fibre connector on the laser console.
- Activate SERVICE MODE 1 using the service plug or the appropriate DIL switch.
- Select Power Mode (PM).
- Go to "Power" and set a power output of 1 W.
- Set a pulse duration of 999 (cw).
- Switch to **Treatment/READY** mode.
- Operate the foot switch.
- Record the value P shown on the power measurement system and the values REAL & MON displayed on the laser display.
- If $P = 1 \text{ W}$,
 $REAL = 1 \pm 0.1 \text{ W}$
 $MON = 1 \pm 0.1 \text{ W}$
 then the calibration is OK. If not, proceed as follows:
 - Release the foot switch.
 - Remove the laser console cover.
 - Remove the cover from the beam deflection unit.

Warning

Unprotected radiation exposure, wear laser safety goggles!

- Operate the foot switch.
- Calibrate the laser using the two step potentiometers (coarse adjustment in 5 steps) and the two helical potentiometers (fine adjustment).
- Adjust potentiometer REAL_GN until the power measurement system displays a power of $P = 1 \pm 0.1 \text{ W}$.
- Adjust the other potentiometers for MON_GN until the laser display displays the same value for MON_GN and REAL_GN.
- Release the foot switch.
- Remove the measuring head with the adapter and the dummy plug.
- Reattach the beam deflection unit and laser console covers.

Calibrating the applicators

- Connect the applicator to be calibrated.
- Switch to SERVICE MODE 1 (SM1).
- Activate the aiming beam. Place the measuring head of the power measurement system behind the applicator so that about 1/3 to 2/3 of the detector surface is illuminated.
- Switch to SERVICE MODE 2.
- Select the currently connected applicator.
- Press the **STANDBY/READY** key to switch to **Treatment/READY** mode.
- Press the foot switch to fire the laser and read off the power.
- If the delivered power downstream of the applicator is more than 1 W, release the foot switch. Increase the applicator transmission factor by 1%. Press the foot switch again and read off the power. Repeat this procedure until the power downstream of the applicator is $1 \text{ W} \pm 2 \%$.
- If the delivered power downstream of the applicator is less than 1 W, release the foot switch. Reduce the applicator transmission factor by 1%. Press the foot switch again and read off the power. Repeat this procedure until the power downstream of the applicator is $1 \text{ W} \pm 2 \%$.
- Return to the main menu and disconnect the service plug.

Ordering information, accessories and spare parts

Component overview:	Catalogue No.
VISULAS 532s	000000-1149-623
LSL 532s laser slit lamp	000000-1149-675
LSL YAG III Combi laser slit lamp	000000-1172-069
Parallel tube f=140 mm	319770-9011-000
Eyepiece 10x	319770-9110-000
VISULINK 532/U (2.0 m)	000000-1137-688
VISULINK 532/U (3.5 m)	000000-1137-689
532-Adapter Haag-Streit 900 [®] BM, 900 [®] BQ	000000-1137-681
532-Adapter Zeiss SL 120 / SL 130	000000-1137-682
532-Adapter Zeiss 20 SL / 30 SL	000000-1137-683
532-Adapter Zeiss SL 115 Classic	000000-1137-684

Endo probes

Sterile endo probes, 20 gauge, straight (5 pcs.)	000000-1212-036
Sterile endo probes, 20 gauge, angled (5 pcs.)	000000-1212-331

Physician's safety filter for surgical microscope

Simple physician's safety filter, active, for Zeiss microscopes	000000-1336-953
Double physician's safety filter, active, for Zeiss microscopes	000000-1338-594
Simple physician's safety filter, fixed, for Zeiss microscopes	000000-1272-301
Simple physician's safety filter, active, for Leica/Wild microscopes	000000-1338-593
Double physician's safety filter, active, for Leica/Wild microscopes	000000-1338-595
Simple physician's safety filter, fixed, for Leica/Wild microscopes	000000-1272-417

Indirect ophthalmoscope

LIO 532	000000-1041-845
Cold light source (230 V)	300746-8003-000
Cold light source (115 V)	300746-8003-710
Halogen lamp 24 V / 150 W	380079-9260-000
Optical fibre	000000-1177-999
Illumination optical fibre	300746-8002-000

Aspherical lenses

The LIO 532 laser indirect ophthalmoscope can be used with all standard aspherical lenses.

To obtain the best results, we recommend the following special aspherical lens, which has a coating optimised for argon and diode laser radiation.

20 D	000000-0400-625
------	-----------------

Contact lenses

Mainster Standard 90° contact lens	306877-9005-000
Mainster Wide Field 125° contact lens	000000-0233-961
Three mirror contact lens	306877-9006-000
Three mirror contact lens, small	000000-0235-830
Ritch trabeculoplasty lens	306877-9020-000

Accessories

Laser safety goggles Argon / 532s in accordance with DIN EN 207	000000-0450-051
Laser safety goggles YAG in accordance with DIN EN 207, IR 1064 L5	000000-1268-052
Transport case (VISULAS 532s with VISULUNK 532/U)	000000-1172-722
Halogen lamp 12 V, 30 W for laser slit lamp	000000-0120-704
Table mounting for VISULAS 532s on IT 4L	000000-1296-221
Table mounting for VISULAS YAG III Combi on IT 4L	000000-1268-416

Carl Zeiss Meditec and their authorised representatives will be happy to provide you with information on further accessories:

- ☐ DigiCam adapter
- ☐ Video documentation equipment
- ☐ Applanation tonometer

VISULAS 532s laser system

Laser type	Frequency doubled Nd:YVO ₄ , diode-pumped, cw	
Laser class	Class 4 (EN 60825-1:2003)	
Wavelength of therapy beam	532 nm	
Pulse duration	10 ... 2500 ms, cw	
Max. power at cornea	1500 mW (slit lamp), 1500 mW (endoprobe)	
Max. power at resonator	3 W	
Pulse interval	Continuously adjustable from 100 to 6000 ms	
Aiming beam	Diode, 620 to 650 nm, max. 1 mW at the cornea, max. 5 mW at the diode	
Environmental conditions	Ambient temperature:	10 to 35 °C
	Relative humidity:	0 to 70 % (non condensing)
	Atmospheric pressure:	700 to 1060 hPa
Storage and shipment conditions	Temperature:	-25 to 55 °C
	Relative humidity:	0 ... 60 %
	Atmospheric pressure:	700 to 1060 hPa
Rated voltage	100 V to 240 V AC (90 V to 264 V)	
Rated frequency	50 / 60 Hz	
Max. current	7 A	
Power consumption	400 VA	
Protection type	IP 20	
Protection Class	I	
Device type	B (in accordance with DIN EN 60601-1)	
Earth conductor	This instrument should only be connected to properly earthed power sockets.	
Cooling	Thermoelectric	
Laser console dimensions	H x W x D = 151 mm x 289 mm x 400 mm	
Laser console weight	15 kg (incl. control panel)	
NOHD	35 m (with LSL) 10 m (with Endo probe) 60 m (with LIO 532)	

LSL 532s laser slit lamp

Laser treatment spot size	Continuously adjustable from 50 to 1000 μm (without contact lens) parfo- cal, larger spot sizes depending on contact lens used
Laser beam delivery	Interlaced with slit illumination system
Illumination	12 V, 30 W brightness continuously adjustable.
Slit adjustment	Slit length variable in steps of: 1 / 3 / 5 / 9 / 14 mm; Slit width continuously adjustable: 0 to 14 mm Slit image rotation: 0°, $\pm 45^\circ$, 90°.
Magnification	5 magnifications, in steps of 5x, 8x, 12x, 20x, 32x
Electrical power supply	Electrical power is supplied via the VISULAS 532s laser console.
Protection type	IP 20
Weight	11 kg
Accessories	Tonometer, assistant's scope, video documentation equipment etc. from the range of accessories for the SL 115, SL 120 and SL 130 slit lamps

LSL YAG III Combi laser slit lamp

Dimensions (incl. laser head)	H x W x D = 623 mm x 350 mm x 400 mm
Electrical power supply	Electrical power is supplied via the VISULAS YAG III laser console.
Weight	12 kg

VISULINK 532/U

Laser treatment spot size	Continuously adjustable from 50 to 1000 μm (without contact lens) parfo- cal, larger spot sizes depending on contact lens used
Dimensions	H x W x D = 180 mm x 35 mm x 125 mm
Weight	0.4 kg

LIO 532 indirect ophthalmoscope

Laser spot size on the fundus of a pati- ent with normal vision using a 20 D asphe- rical lens	360 μm (Aiming and therapy beam)	
Laser optical fibre	Length	300 cm
	Core diameter	160 μm
Illuminating optical fibre	Length	300 cm
	Fibre bundle diameter	4 mm
Physician's safety filter	Transmission < 0.005 % for 532 nm	

Technical details and delivery package subject to change.

This instrument meets the requirements of the EU Medical Device Directive 93/42/EEC and its national implementation in the form of the German Medical Products Act (MPA).

Device class under MPA: II b

UMDNS-No.: 17-446 (VISULAS 532s)

UMDNS-No.: 12-818 (LIO 532)

UMDNS-No.: 12-281 (Slit lamp)

This declaration is rendered invalid if changes are made to the product without the manufacturer's authorisation.

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Carl Zeiss Meditec AG