Linvatec[®], Hall[®] Surgical, Concept[®] and Shutt[®] Manual Instruments

ENGLISH	1
FRANÇAIS	2
ITALIANO	3
DEUTSCH	4
ESPAÑOL	5
ΕΛΛΗΝΙΚΑ	6
PORTUGUÊS	7
PORTUGUÊS DO BRASIL	8
NEDERLANDS	9
DANSK	10
SVENSKA	11
SUOMI	12
NORSK	13
简体中文	14
日本語	15
POLSKI	16
TÜRKÇE	17
ROMÂNĂ	18
العربية	19
РУССКИЙ	20
HRVATSKI	21
ČEŠTINA	22
MAGYAR	23
LIETUVIŲ K.	24
SLOVENČINA	25
SLOVENŠČINA	26
한국어	27

Linvatec®, Hall® Surgical, Concept® and Shutt® Manual Instruments

CLEANING AND STERILIZATION GUIDELINES

NOTE: Prior to use, cleaning and sterilization, remove all protective packaging and tip protector.

GENERAL

Inspect instruments before and after processing for proper function and alignment of pins and for chipping of plated surfaces.

CLEANING

Open and/or disassemble instrument immediately after use. Place in a basin with distilled water.

NOTE: Saline solution is not recommended as it has a corrosive effect on certain metals.

Instruments should be cleaned using a neutral pH cleaning solution and a non-abrasive brush.

A commercial water-based lubricant should be used to minimize wear of movable parts.

STEAM STERILIZATION

Prepare instruments such that all surfaces have direct contact with steam. Hinged instruments and luer-locks should be in the open position. Disassemble the instruments whenever possible. Double wrap instruments using FDA cleared sterilization wrap.

Method	Cycle	Minimum Temperature	Minimum Exposure	Minimum Dry Cycle
Steam (Double Wrapped)	Pre-vacuum	270°F (132°C)	3 minutes	8 minutes
Steam (Double Wrapped)	Gravity	250°F (121°C)	30 minutes	8 minutes

For sterilization of instruments in a tray, consult guidelines provided with the specific sterilization tray.

NOTE: Guidelines provided above do not guarantee sterility. Your institution is responsible for the normal sterility assurance validation.

Reference: AAMI (Association for the Advancement of Medical Instrumentation.)

CONTACT INFORMATION

For more information or a product demonstration, contact your regional ConMed Linvatec sales representative, or call 1-866-426-6633 in the U.S.

Instruments manuels Linvatec®, Hall® Surgical, Concept® et Shutt®

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All rights reserved. P000017745 Rev B 03/2022
MDSS GmbH
CONMED UK Ltd