

---

**THE BOC GROUP**

# **Biox 3760 Pulse Oximeter**

---

## **Operation and Maintenance Manual**

## **Important**

This manual is subject to periodic review, update, and revision. Customers are cautioned to verify that the manual's information applies to the software and hardware present in the equipment.

This product performs as described in this manual, and in accompanying labels and/or inserts, when assembled, operated, maintained, and repaired in accordance with the instructions provided.

This product must be cleaned and checked periodically. Do not use a defective product. Parts that are broken, missing, plainly worn, distorted, or contaminated should be replaced immediately. If repair or replacement become necessary, call or write to request service advice from the nearest Ohmeda Regional Service Center (listed on back cover). Do not repair this product or any of its parts other than in accordance with written instructions provided by Ohmeda and by Ohmeda-trained personnel.

The product must not be altered without the prior written approval of Ohmeda's Safety Department. The user of this product shall have the sole responsibility for any malfunction that results from improper use, faulty maintenance, improper repair, unauthorized service, damage, or alteration by anyone other than Ohmeda.

The safety, reliability, and performance of this device can only be assured under the following conditions:

- If the device has been used according to the accompanying operating instructions.
- If fittings, extensions, readjustments, changes, or repairs have been carried out by Ohmeda's authorized agents.
- If it is used in buildings having ground equalization wiring that complies with relevant IEC or local standards and regulations (UL, ETL, CSA, BSI, TUV, etc.).

**Ohmeda is a registered trademark of The BOC Group.**

All specifications subject to change without notice.

Text revised May 1994

<b>Table of Contents</b>	i
<b>List of Figures</b>	iii
<b>1/Overview</b>	
Introduction	1-1
Description	1-1
Principles of operation	1-1
Theory	1-1
Functional components	1-3
Features and controls	1-4
Front Panel	1-4
Rear Panel	1-5
Precautions	1-5
Warnings	1-5
Cautions	1-7
<b>2/Operation</b>	
General operation guidelines	2-1
Preoperative checklist	2-1
Visual inspection	2-1
Functional inspection	2-2
Data validity	2-2
Clock set mode	2-3
Printer operation	2-4
Side panel—printer door	2-4
The print-out	2-7
Switching formats/speeds while printing	2-8
Interrupting printing	2-9
Stopping printing	2-9
Statistics	2-9
Loading printer paper	2-10
<b>3/Messages and Troubleshooting</b>	
Low SpO <sub>2</sub> alarm limit violation	3-1
Probe alarm conditions	3-1
Probe alarm messages	3-2
Signal and device failure problems	3-3
<b>4/Maintenance and Service</b>	
Cleaning	4-1
Oximeter	4-1
Probes	4-2
Recharging the battery	4-2
Repair policy and procedure	4-2
Obtaining service	4-3
Packaging and return procedure	4-3
Accessories	4-4
<b>A/Specifications</b>	A-1
References	A-4
<b>Index</b>	Index-1
<b>Warranty</b>	W-1

<b>Figure No.</b>	<b>Description</b>	<b>Page</b>
1-1	Comparative light absorbtion.....	1-2
1-2	Extinction vs. wavelength graph .....	1-2
1-3	Functional components .....	1-3
1-4	Front Panel.....	1-4
1-5	Rear Panel .....	1-5
2-1	Side panel—printer door.....	2-4
2-2	Numeric format.....	2-5
2-3	Graphic format.....	2-6
2-4	Graphic resolution print-out.....	2-6
2-5	Graphic format (filled and open).....	2-7
2-6	Print-out.....	2-7
2-7	Print-out (continued) .....	2-8
2-8	Loading printer paper .....	2-10

# 1/Overview

## Introduction

This manual describes the proper operation and maintenance for the Ohmeda Biox 3760 Pulse Oximeter (software Revision B). Please read this manual before using the pulse oximeter, paying attention to all details of correct operation along with precautionary measures recommended. You may perform the maintenance procedures included in this manual.

This chapter contains

- A general description of the pulse oximeter, including the theory of operation and functional components.
- An overview of the unit's features and controls.
- A list of all precautions you must follow when using the pulse oximeter.

## Description

The Ohmeda Biox 3760 Pulse Oximeter is a stand-alone, non-invasive, arterial oxygen saturation monitor. It provides continuous, real-time SpO<sub>2</sub> and pulse rate readings.

## Principles of operation

### Theory

The 3760 pulse oximeter operates on the basic assumption that hemoglobin exists in the blood in two principal forms :

- **Oxygenated (HbO<sub>2</sub>)**—with loosely bound O<sub>2</sub> molecules.
- **Reduced (Hb)**—with no bound O<sub>2</sub> molecules.

Arterial oxygen saturation (SpO<sub>2</sub>) is defined as the ratio of oxygenated hemoglobin (HbO<sub>2</sub>) to total hemoglobin (HbO<sub>2</sub> + Hb):

$$SpO_2 = \frac{HbO_2}{HbO_2 + Hb + \text{others}^*}$$

\* Others = carboxyhemoglobin, methemoglobin, sulfhemoglobin, etc.; also see information about interfering substances in the A/Specifications.

**Note:** Because the presence of dyshemoglobins or other pigments cannot be measured by a two-wavelength instrument, the standard, internationally recognized definition for the  $\text{SaO}_2$  read by pulse oximeters is  $\text{SpO}_2$ . The presence of appreciable amounts of these substances may result in erroneous readings. The information in this manual applies whether your unit and the messages it generates refer to  $\text{SaO}_2$  or  $\text{SpO}_2$ .

The attenuation of light energy due to arterial blood flow can be detected and isolated. The oximeter, relying on the observation that arterial blood flow pulsates and other fluids and tissues do not, uses a patented two-wavelength, pulsatile system to differentiate between the light absorbed by hemoglobin and that absorbed by other fluid and tissue constituents. The arterial blood flow pulsation modulates light passing through the blood; nonpulsing fluids and tissues do not modulate the light.




Figure 1-1. Comparative light absorption





Figure 1-2. Extinction vs. wavelength graph

As shown above,  $\text{HbO}_2$  and  $\text{Hb}$  exhibit markedly different absorption (extinction) characteristics to red light at 660 nm and infrared light at 940 nm. The oximeter measures the relative absorption by  $\text{HbO}_2$  and  $\text{Hb}$  at these two wavelengths, and then converts this relative light intensity information into  $\text{SpO}_2$  and pulse rate values.

## Functional components

The Ohmeda Biox 3760 Pulse Oximeter uses electrical components to determine SpO<sub>2</sub> and pulse rate values. The key elements are:

- the probe
- the processing of the probe signal
- the calculations made by the microprocessor



**Figure 1-3.** Functional components

The probe consists of:

- the light source - a red LED and an infrared LED
- the photodetector - an electronic device that produces an electrical current proportional to incident light intensity

The two wavelengths of light generated by the LEDs pass through the tissue at the probe site. This light, which is partially absorbed and modulated, is then collected by the photodetector and converted into an electronic signal. This signal is sent to the oximeter for further processing.

The electronic circuitry takes the current generated by the photodetector, processes it, and passes it to the micro-processor for calculation of the SpO<sub>2</sub> and pulse rate.

The calculation of SpO<sub>2</sub> assumes 1.6% carboxyhemoglobin, 0.4% methemoglobin, and no other pigments. Appreciable variation from these values will influence the accuracy of SpO<sub>2</sub>. These values are based on the Ohmeda Biox 3700 Empirical Calibration Study.

The microprocessor calculates the SpO<sub>2</sub> 30/25 (60/50 Hz) times per second. The calculations are averaged by a running weighted average\* method to determine the displayed SpO<sub>2</sub>. The displayed average uses data over a 6 second period and is updated every 0.67 seconds (60 Hz) or 0.5 seconds (50 Hz).

\* Obtained by assigning a weight (value) to each calculation based on the signal strength and the current average saturation.

The running weighted average method allows erroneous SpO<sub>2</sub> values to be discarded from the determination of the displayed SpO<sub>2</sub>. Erroneous values result from probe movement, electrosurgery, and other sources of signal interference. This method of averaging provides a stable reading, with low sensitivity to interference while retaining the capability of responding quickly to saturation changes.

## Features and controls

### Front panel




Figure 1-4. Front panel

#### 1 Power key

Pressing this turns the oximeter ON or OFF.

#### 2 Alarm Limit key

Adjusts low SpO<sub>2</sub> alarm limit by pressing the key until desired limit is displayed and then release. (Default = 90%)

#### 3 Alarm Silence key

Silences all audible alarms for 60 seconds, regardless if other alarms occur during this 60 second interval.

When pressed, the flashing red alarm light changes to a steady red light to indicate alarm silence. If an alarm condition still exists after the 60 second alarm silence period, the audible tone and flashing light resume.

**Exception:** For OFF PATIENT or NO PROBE alarms, the Alarm Silence key silences the audible alarm until either:

- The specific alarm condition is remedied or
- A different alarm condition is detected.

#### 4 Liquid crystal display

Displays SpO<sub>2</sub>, pulse rate, status messages, and alarm messages.

#### 5 Probe plug connector

The probes supplied with this model of oximeter plug into this nine-hole connector.



**X Ball rod**

Folds half-way out from beneath the oximeter chassis to raise the position of the display for better viewing.

**Important:** Do not fold the rod FORWARD past the half-way stop, as damage to the chassis will occur.


**Rear panel****6 Charger connector**

Figure 1-5. Rear panel

**6 Charger connector**

Connects the oximeter to the AC charger for continuous operation and/or charging of the battery.

## Precautions

Two types of precautions appear in this manual: warnings and cautions.

**Warnings**

A **WARNING** indicates the possibility of injury to the patient or operator.

**Handle the monitor with care**

Improper handling can cause damage or inaccurate results.

**Failure of operation**

If the oximeter fails any part of the preoperative checkout procedures, calibration, or current leakage test, remove it from operation until qualified service personnel have corrected the situation.

The oximeter is a microprocessor-based device designed to immediately shut down if the microprocessor fails. This prevents the possible display of erroneous information. No alarms forewarn this action.

**Data validity**

Calibration is verified during powerup. Do not operate the oximeter unless it is properly calibrated or inaccurate patient readings will result.

Excessive ambient light, excessive motion, low perfusion, or electrical interference at the probe site may cause erroneous readings.

To prevent erroneous readings, do not use an inflated blood pressure cuff on the same limb as the oximeter probe.

### **Electrical shock hazard**

This equipment must be properly grounded.

- Connect this equipment only to a three-wire, grounded, hospital-grade receptacle. The three-connector plug must be inserted into a properly wired three-wire receptacle; if a three-wire receptacle is not available, a qualified electrician must install one in accordance with the governing electrical code.
- Do not under any circumstances remove the grounding connector from the power plug.
- Do not use extension cords or adapters of any type. The power cord and plug must be intact and undamaged.
- If there is any doubt about the integrity of the protective earth conductor arrangement, operate the monitor on internal battery power until the AC power supply protective conductor is fully operational.

Do not remove the monitor cover. An operator may only perform maintenance procedures specifically described in this manual. Refer servicing to qualified service personnel trained in the repair of this equipment.

Use **only** the charger supplied for this oximeter.

### **Electrical shock and flammability hazard**

Before cleaning the oximeter, always turn it off and disconnect the power cord from the AC power supply.

### **Explosion hazard**

Do not use the oximeter in the presence of flammable anesthetics or other flammable substances.

### **Battery replacement**

Unauthorized personnel should not attempt to install, connect, or replace the oximeter's battery.

- Removing the cover and/or faulty battery connections could be hazardous and will void the warranty.
- Reversing the battery connections could result in injury and will permanently damage the circuitry.
- If trained technical personnel are not available, call Ohmeda for assistance.
- For proper operation, replace only with the appropriate Ohmeda battery.

### **Patient safety**

To prevent patient injury or equipment damage, use only oximeter probes identified for use with this monitor (see the instructions for the probe you are using).

If a probe is damaged in any way, discontinue use immediately.

Prolonged monitoring or patient condition may require periodically changing the probe test site. To reduce the risk of blistering, skin erosion, or ischemic skin necrosis, change the probe site as specified in the user instructions for the probe or sensor you are using. If any evidence of the above conditions appears before the specified time (for example, discoloration or reddening), change the probe site immediately.

To avoid any possibility of patient discomfort or injury during magnetic resonance imaging (MRI),

- Do not allow the oximeter probe cable to come in contact with the patient's body; keep the cable off of the patient or place a blanket or other insulating material between the patient and the probe cable.
- Position the oximeter probe and probe cable as far from the center of the magnetic field as possible.

The correct use of the oximeter is to measure only arterial oxygen saturation (SpO<sub>2</sub>) and pulse rate.

- A pulse oximeter does not measure respiration and under no circumstances should be used as a substitute for an apnea monitor.
- The oximeter must not be used as the primary monitor for infants being monitored for apnea, either in the hospital or in the home setting. It measure SpO<sub>2</sub> and pulse rate, and only in conjunction with other appropriate monitoring techniques.
- A pulse oximeter is often used during sleep studies with adults, but must be used only to gather information regarding SpO<sub>2</sub> and pulse rate during these studies.
- A pulse oximeter is to be used only by or on the order of medically trained personnel.

Refer to the instructions for the probe you are using for detailed warning information with respect to probes.

## Cautions

A **CAUTION** indicates a condition that may lead to equipment damage or malfunction.

Federal law in the USA and Canada restricts the sale of this device by or on the order of a licensed medical practitioner.

Always make sure the monitor is set up to operate at the AC power supply voltage present at the "wall" receptacle.

Avoid storing the oximeter and probes outside the following temperature ranges:

- Units without paper: -20° C (-4° F) to 60° C (140° F)
- Units with paper: -20° C (-4° F) to 45° C (115° F)

To prevent damage to the lead-acid battery, do not turn the monitor on after the **Recharge Battery** alert message appears without first plugging it into the AC power supply.

Do not operate the printer without paper. Damage to the print head may result. Paper purchased from Ohmeda has red edges to indicate the end of the paper on the roll is approaching.

### **Cleaning**

- Do not autoclave, pressure sterilize, or gas sterilize this oximeter.
- Do not soak or immerse the monitor in any liquid.
- Use the cleaning solution sparingly. Excessive solution can flow into the monitor and cause damage to internal components.
- Do not touch, press, or rub the display panel with abrasive cleaning compounds, instruments, brushes, rough-surface materials, or bring it into contact with anything that could scratch the panel. Do not use solutions containing acetone or other harsh solvents to clean the display panel.

Detailed information for more extensive repairs is included in the service manual solely for the convenience of users having proper knowledge, tools, and test equipment, and for service representatives trained by Ohmeda.

Refer to the instructions for the probe you are using for detailed cleaning and caution information with respect to probes.

## 2/Operation

### General operation guidelines

This chapter contains

- The preoperative checklist
- Information to help you assure validity of the pulse oximeter's data.
- Instructions for setting the oximeter's internal clock.
- Instructions for operating the printer.

#### **WARNINGS:**

**Data validity—Excessive ambient light, excessive motion, low perfusion, or electrical interference at the probe site can result in the display of invalid data.**

**Electrical shock hazard—Use only the charger supplied for this oximeter. (CAUTION: Damage to the equipment may also result.)**

**Explosion hazard—Do not use in the presence of flammable anesthetics or other flammable substances.**

**Failure of operation—If the oximeter fails to respond as described, do not use it until the situation has been corrected by Ohmeda-trained service personnel.**

**Patient safety—If a probe is damaged in any way, discontinue use immediately.**

#### **CAUTIONS:**

Use **only** Ohmeda probes, sensors, and interconnect cables with this oximeter. Otherwise, damage to the equipment may result.

Do not apply tension to the probe cable. Probe damage may result.

### Preoperative checklist

Perform the following tests daily to assure proper operation of the oximeter and probes.

#### **Visual inspection**

1. Inspect the oximeter case for damage.
2. Ensure the display windows are clean. See 4/Maintenance and Service, Cleaning.

### Functional inspection

1. Plug the charger into AC mains power and into the charging connector on the rear panel of the oximeter.
2. Check the probe.
  - a. If appropriate, make sure the probe opens and closes smoothly. If there is any unevenness or variation in the closing motion, replace the probe.
  - b. Ensure that there is no foreign material such as tape or cotton covering either the emitter or detector.
3. Connect an Ohmeda probe to the oximeter. Follow the application instructions carefully for the probe or sensor you are using.

Check that the probe connector makes a firm connection with the oximeter.

Check that the probe cable is not twisted, sliced, or frayed.

4. Attach the probe to either a finger or an ear.
5. To turn the oximeter ON, press the **Power** key. Verify that the oximeter displays:

\* POWER ON \*

then runs the self-diagnostic test and displays:

SYSTEM CHECK ...

In a few seconds, SpO<sub>2</sub> and pulse rate data appear.

**Important:** If the oximeter fails any part of the diagnostic self-test during power up, the message SERVICE UNIT appears. Turn the oximeter off and refer to Ohmeda-trained service personnel.

6. Verify that the alarms are functional:
  - Remove the probe from the finger or ear. Ensure that the OFF PATIENT alarm message appears, the alarm tone sounds, and the red alarm light flashes.

**Note:** A different OFF PATIENT message may appear when using the Flex II Probe or the SoftProbe.
  - Unplug the probe from the oximeter. Ensure that the NO PROBE alarm message appears, the alarm tone sounds, and the red alarm light flashes.
7. Turn the oximeter off. The display blanks.

### Data validity

The probe must be correctly attached to the patient and the data must be reliable.

**WARNING:** Data validity—Excessive ambient light, excessive motion, low perfusion, or electrical interference at the probe site can result in the display of invalid data.

Factors affecting reliability of the signal and possible display of inaccurate SpO<sub>2</sub> readings are:

- Motion at the probe site
- Electronic noise
- Signal strength

These may be detected when the oximeter's pulse rate varies significantly from the palpated pulse rate. If you suspect that this may be the case, check the following:

- **Excessive patient movement or detector not flush with the skin:** Be sure that the probe is properly attached. If necessary, tape the probe and/or probe cable.
- **Electronic noise:** The electronic noise emitted from electrosurgery, MRI (NMR), and some electrical devices such as fans can interfere with the signal from the probe.

**WARNING: Patient Safety—To avoid any possibility of patient discomfort or injury during magnetic resonance imaging,**

- **Do not allow the oximeter probe cable to come in contact with the patient's body; keep the cable off of the patient or place a blanket or other insulating material between the patient and the probe cable.**
- **Position the oximeter probe and probe cable as far from the center of the magnetic field as possible.**
- **Probe site too thick:** A large distance between the emitter (LED) and the detector can reduce the signal strength and result in a poor signal. Select a thinner site.
- **Poorly perfused site:** Select an alternate site or rub the site to increase perfusion.
- **Dark pigmentation:** Very dark pigmentation may result in low signal strength. Move the probe to a site with less pigmentation (e.g., the finger).
- **Artificial fingernails or heavy fingernail polish:** Remove polish or artificial fingernail, or select another probe site.

### Clock set mode

1. Turn the oximeter on. To be sure the date and time are accurately printed, the internal clock must be set correctly.
2. Hold down the **Stop** key (for approximately 3 seconds) until the front panel displays the date and time in the following format:

mm/dd/yy hh:mm

3. To change any of the clock settings, press the key that corresponds to the section you wish to change as shown in the following chart:

Key:	Print Slow	Print Fast	Graph Slow	Graph Fast	Pause Resume
Sets:	Month	Day	Year	Hour	Minute

- Pressing a key once changes the setting by one increment.
  - Pressing a key continuously changes increments until you remove your finger from the key.
4. When you are done setting the clock, press the **Stop** key. The oximeter resumes its normal operation.

**Important:** If an alarm condition occurs while you are in the clock set mode, the audible alarm sounds and the red alarm light flashes, but the message will not appear on the display. (Should you want to see the alarm condition message, press the **Stop** key to exit the clock set mode.)

### Printer operation

If you are using an Ohmeda Biox 3760 Pulse Oximeter with a built-in printer, read this section.

**CAUTION:** Do **not** operate the printer without paper. Damage to the print head may result. (Paper purchased from Ohmeda will have red edges to indicate the end of the roll.)

- The printer records the patient's SpO<sub>2</sub> and pulse rate readings in one of six different formats (user-selectable).
- If an alarm condition occurs, the corresponding message is printed. (If a LO BAT condition occurs, the printer automatically records the closing statistics and shuts down.)
- Thermal paper for the printer is available from Ohmeda (see "Accessories" in 4/Maintenance and Service). Use of other manufacturers' paper is **not** recommended as damage to the printer mechanism may occur.
- If long-term (more than 2 years) storage of a print-out is desired, a photocopy should be made and stored instead of the original. (Thermal paper tends to fade after 2-3 years.)
- Storage temperature range (thermal paper): -20° C (-4° F) to 45° C (115° F).

#### Side Panel—Printer Door




Figure 2-1. Side panel—printer door



The Ohmeda Biox 3760 Pulse Oximeter with Printer records patient SpO<sub>2</sub> and pulse rate in either a numeric format or graphic format.

1 through 4: Press the key that selects the printed output format that you want, as described below.

### Numeric format

#### Print slow

Prints the lowest SpO<sub>2</sub> and its corresponding pulse rate reading over a 30-second period, every 30 seconds.

(For 50 Hz Models, substitute 24 seconds for 30.)

#### Print fast

Prints the lowest SpO<sub>2</sub> and its corresponding pulse rate reading over a 6-second period, every 6 seconds.

(For 50 Hz Models, substitute 4 seconds for 6.)

```

NAME: .....
COMMENTS: .....
.....
.....
BY: .....
DATE - SEP/17/93
TIME - 10:44:07
PRINT FAST, 6 SEC./READING

SaO2%  PR    TIME
94%    61    10:44:13
94%    60    10:44:19
94%    60    10:44:25
94%    60    10:44:31
94%    60    10:44:37
94%    61    10:44:43
93%    65    10:44:49
88%    79    10:44:55

```

Figure 2-2. Numeric format

### Graphic format

#### Graphic slow

Prints data grouped in 48-second periods (broken into eight 6-second intervals). Each line of dots represents the lowest SpO<sub>2</sub> per interval.

The time each 48-second period begins is printed to the left of the graph. The pulse rate measured at the time of the lowest SpO<sub>2</sub> for that same period is printed to the right of the graph.

(For 50 Hz Models, substitute 32 seconds for 48 and 4 seconds for 6.)

```

NAME: .....
COMMENTS: .....
.....
.....
BY: .....
DATE - SEP/17/93
TIME - 10:50:41
GRAPH SLOW, 6 SECONDS/DOT
HR:MIN      SaO2%      PR
50 60 70 80 90 100
10:50 [graph] 81
10:51 [graph] 71
10:52 [graph] 74
10:53 [graph] 79
10:53 [graph] 82

```

Figure 2-3. Graphic format

**Graphic fast**

Prints data grouped in 16-second periods (broken into eight 2-second intervals). Each line of dots represents the lowest SpO<sub>2</sub> per interval.

The time each 16-second period begins is printed to the left of the graph. The pulse rate measured at the time of the lowest SpO<sub>2</sub> for that same period is printed to the right of the graph.

(For 50 Hz Models, substitute 10.4 seconds for 16 and 1.3 seconds for 2.)

```

NAME: .....
COMMENTS: .....
.....
.....
BY: .....
DATE - SEP/17/93
TIME - 11:00:12
GRAPH SLOW, 6 SECONDS/DOT
HR:MIN      SaO2%      PR
50 60 70 80 90 100
11:00 [graph] 83
11:01 [graph] 74

```

Figure 2-4. Graphic resolution print-out

### Filled graph or open graph

The printer automatically prints a filled graph when either the **Graph Slow** or the **Graph Fast** key is pressed.




Figure 2-5. Graphic format (filled and open)

To change between a filled graph and an open graph:

1. Hold down the **Graph Fast** key (for approximately 3 seconds) until the front panel displays the status message **OPEN GRAPH** or **FILLED GRAPH**.
2. To clear that status message and return to the present patient  $\text{SpO}_2$  and pulse rate reading, press the **Stop** key.
3. To start printing, press the desired **Graph** key.

### The print-out




Figure 2-6. Print-out

- 1 Date and time the printer was turned on (set from the oximeter's internal clock).
- 2 Print speed and format selected, and the time frame per reading.
- 3 SpO<sub>2</sub> and corresponding pulse rate and time (shown as SpO<sub>2</sub> and in the numeric format).
- 4 Alarm messages that print as the condition occurs.

## 5 → \*\* 20 MINUTE STATISTICS \*\*

```

      CUMULATIVE
      0% OF TIME SaO2 < 90%
      0% OF TIME SaO2 < 85%
      0% OF TIME SaO2 < 80%
      0% OF TIME SaO2 < 70%
      LOWEST SaO2 WAS -
      97% AT 11:33:25
      AVERAGE SaO2 = 98%

```

```

      BY RANGE
      100% OF TIME SaO2 90 - 100%
      0% OF TIME SaO2 85 - 89%
      0% OF TIME SaO2 80 - 84%
      0% OF TIME SaO2 70 - 79%
      0% OF TIME SaO2 < 70%

```

## HISTOGRAM

```

      % TIME PER SaO2 RANGE
      0 20 40 60 80 100 %
      90-100 ██████████ 100
      85- 89 ██████████ 0
      80- 84 ██████████ 0
      70- 79 ██████████ 0
      < 70 ██████████ 0

```

## PRINT FAST, 6 SEC./READING

```

      SaO2% PR TIME
      97% 76 11:54:01
      97% 76 11:54:07
      97% 76 11:54:13
      97% 76 11:54:19
      97% 76 11:54:25

```

## 6 → \*\*\* CLOSING STATISTICS \*\*\*

```

      CUMULATIVE
      0% OF TIME SaO2 < 90%
      0% OF TIME SaO2 < 85%
      0% OF TIME SaO2 < 80%
      0% OF TIME SaO2 < 70%
      LOWEST SaO2 WAS -
      97% AT 11:06:53
      AVERAGE SaO2 = 97%

```

```

      BY RANGE
      100% OF TIME SaO2 90 - 100%
      0% OF TIME SaO2 85 - 89%
      0% OF TIME SaO2 80 - 84%
      0% OF TIME SaO2 70 - 79%
      0% OF TIME SaO2 < 70%

```

## HISTOGRAM

```

      % TIME PER SaO2 RANGE
      0 20 40 60 80 100 %
      90-100 ██████████ 100
      85- 89 ██████████ 0
      80- 84 ██████████ 0
      70- 79 ██████████ 0
      < 70 ██████████ 0

```


**OhmedaBiox**   
**3760 Pulse Oximeter**

Figure 2-7. Print-out (continued)

- 5 Trend statistics that automatically print out during monitoring.
  - Every 20 minutes (print fast or graph fast)
  - Every 60 minutes (print slow or graph slow).
- 6 Closing statistics that print when
  - the Stop key or another print key is pressed.
  - the LO BAT condition occurs.
  - the oximeter memory is filled.

**Switching formats/speeds while printing**

To select a different print format or print speed while printing, press the key corresponding to the selection. The oximeter prints the closing statistics and starts printing in the newly selected print format or speed.

### Interrupting printing

1. To make the printer pause, press the **Pause/Resume** key once. The printer finishes printing the current line (if printing) and then prints:

----- PAUSE -----

The SpO<sub>2</sub> and pulse rate values continue to be displayed on the front panel (provided a probe is connected to the oximeter and patient). The SpO<sub>2</sub> and pulse rate values are not collected during this time for the closing statistics.

2. To resume printing, press the **Pause/Resume** key again. The printer prints:

----- RESUME -----

and then prints in the previously selected format.

### Stopping printing

To stop the printer, press the **Stop** key once. Closing statistics are printed and then the printer turns off. If you are still monitoring a patient, SpO<sub>2</sub> and pulse rate values continue to be displayed on the front panel.

### Statistics

#### Closing statistics

Closing statistics are calculated from the patient SpO<sub>2</sub> readings taken while the printer was activated. They do not include patient SpO<sub>2</sub> readings taken when the pause/resume function was engaged. Closing statistics include:

- Cumulative—contains the percentage of time the patient's SpO<sub>2</sub> was below 90%, 85%, 80%, and 70%. It also contains the time and reading of the lowest SpO<sub>2</sub>; and prints out the average SpO<sub>2</sub> reading.
- By range—shows the percentage of time the SpO<sub>2</sub> was in each of five ranges (90-100%, 85-89%, 80-84%, 70-79%, and < 70%). Due to the way the software calculates and rounds data, it is possible that the total of these percentages may be less than 100%.

(Note: "By range" statistics are only available in Revision B software.)

- Histogram—presents the "by range" statistics in a graphic form and is always printed as a filled graph.

The closing statistics are printed when

- The **Stop** key is pressed.
- Another print format is selected.
- The alarm message LO BAT occurs.
- The oximeter's memory is filled.

### Trend statistics


Trend statistics are calculated the same as closing statistics and contain the same three parts (cumulative, by range, and histogram) as described above.

They are automatically printed out every 20 minutes while in the print fast or graph fast modes, or every 60 minutes while in the print slow or graph slow modes.

### Loading printer paper

**CAUTION:** Never pull the paper backwards through the paper feed, this may damage the feeder mechanism. To remove the paper, cut between the paper roll and the printer door; then pull the paper forward through the front of the printer door.

1. Turn the oximeter on.
2. Place the oximeter so that the printer door can open to its maximum position. To open the printer door, press down on the thumb grip (ribbed indentation) on the top of the printer door and pull back.



**Figure 2-8.** Loading printer paper

3. Unfasten the spindle cap from the paper spindle. Remove the paper roll if present.
4. Slide the new paper roll onto the spindle so that the leading edge of the paper comes out from the beneath the roll.
5. Refasten the spindle cap securely.
6. Make sure the leading end of the paper has a clean, straight edge. Feed the paper into the slot in the printer door, gently guiding it until the feeder mechanism catches it. **Do not force the paper.**
7. Close the printer door.
8. To activate the paper feed, press the **Print Fast** key for approximately 3 seconds until the front panel displays PAPER FEED.
9. When the paper is visible outside of the printer door, press the **Stop** key.

## **3/Messages and Troubleshooting**

This chapter provides information in the event of instrument or signal problems. It is divided into four parts:

- Low SpO<sub>2</sub> alarm limit violation
- Probe alarm conditions
- Probe alarm messages
- Signal and device failure problems, including messages that may appear.

### **Low SpO<sub>2</sub> alarm limit violation**

When the low SpO<sub>2</sub> alarm limit is violated:

- An alarm tone sounds.
- The red alarm light flashes.

Pressing the Alarm Silence key temporarily suppresses the audible alarm for 60 seconds and changes the flashing red Alarm light to steady red.

### **Probe alarm conditions**

These occur when the oximeter detects conditions affecting the probe, its placement, or a probe failure. Alarms can be silenced for 60 seconds with the following exceptions:

- When OFF PATIENT or NO PROBE alarms occur:
  - An alarm tone sounds.
  - The red alarm light flashes.
  - An alarm message appears on the display.
- Pressing the Alarm Silence key during either OFF PATIENT or NO PROBE alarms silences the audible alarm and changes the red alarm light from flashing to steady until
  - The specific alarm condition is remedied.
  - A different alarm condition is detected.
  - A different message is displayed.

#### Probe alarm messages

MESSAGE	POSSIBLE CAUSE(S)	ACTION(S)
ID ERROR	Probe is incompatible with your unit	Use only Ohmeda probes with the correct connector.
LOW LIGHT	Dirt on probe emitter or detector, or at test site.  Tissue sample too thick.  Misaligned or poorly positioned probe.  Emitter light blocked by fingernail polish, artificial fingernails, or dark pigmentation.  Detector failure.	Clean the probe and/or test site.  Select thinner test site.  Reposition the probe or select an alternate test site.  Remove polish or artificial fingernail, or use alternate test site with less coloration.  Replace probe.
LOW QUAL SGNL	Poorly perfused site.  Tissue sample too thick.	Rub site vigorously 20-30 seconds.  Check patient and oximeter setup.  Select thinner site (i.e., less distance between the probe emitter and detector).
NO PROBE	Probe not plugged in or fully inserted into probe connector.	Insert probe plug fully into probe connector.  Use only Ohmeda probes with the correct connector.
OFF PATIENT	Probe off patient.	Attach the probe to the patient.
<b>Note:</b> Other messages may occur with Flex II Probe, SoftProbe, or Easy Probe.	Excessive light detected.  Extremely thin tissue at test site.	Shield probe site from ambient light.  Select an alternate test site.
PROBE/UNIT FAIL  (Oximeter alarms for 2 seconds and then shuts off)	Oximeter detected probe or probe circuit failure.	Replace probe. If the oximeter still won't operate, unit requires service.



## Signal and device failure problems

SYMPTOM/MESSAGE	POSSIBLE CAUSE(S)	ACTION(S)
CHECKSUM ERROR	Printer EPROM error.	Unit requires service.
FIFO OVERFLOW	Buffer overflow.	Turn the oximeter off and back on. If the message persists, unit requires service.
	Paper jammed in printer mechanism (printer motor not running).	Try pulling paper out through the front of the printer door. <b>Do not pull backwards.</b>
INTERFERENCE	Strong radio frequency interference; electro-cautery.	Operate from battery power or plug charger into a different wall outlet (AC mains power).
LO BAT	Battery is getting low.	Recharge the battery or operate from AC mains power. (Note: LO BAT message replaces the pulse rate reading.)
SERVICE UNIT	An internal component has failed.	Turn the unit off; unit requires service.
SYSTEM ERROR X (X = alphanumeric value)	Internal hardware failure.	Note error code and turn oximeter off; unit requires service.
RECHARGE BATT	Battery is unable to supply proper operating voltage.	Recharge the battery or operate from AC mains power. (Note: The unit turns off approximately 10 seconds after RECHARGE BATT message appears.)
Display blank or difficult to read	Display contrast requires adjustment.	Unit requires service.
	Power off.	Press power switch to turn unit on.
	Battery needs recharging or replacement.	Plug oximeter into charger and AC mains power. If the battery won't recharge, unit requires service.
	Blown fuse or internal circuit failure.	Plug oximeter into AC mains power. If the oximeter still will not operate, unit requires service.
Non-correlating or erratic pulse	Excessive patient motion.	Select alternate probe site, or tape probe cable to patient.
Non-correlating SpO <sub>2</sub>	Comparison vs. calculated SpO <sub>2</sub> .	Compare with a co-oximeter.
(i.e., oximeter reading does not correspond with co-oximeter reading)	Excessive ambient light.	Shield probe from ambient light.
	High concentration of carboxyhemoglobin or other dysfunctional hemoglobin.	Use oximeter only to monitor SpO <sub>2</sub> trends, not values.

#### **Notes**

## **4/Maintenance and Service**

This chapter contains the following information:

- How to clean the oximeter and probes
- How to recharge the battery
- How to obtain service and repairs
- A list of accessories you can order for your unit

On a routine basis, you must:

- Clean the exterior surface of the monitor—at least weekly.
- Clean reusable probes—after each patient use.
- Recharge the battery—whenever the LO BAT or RECHARGE BATT message appears and whenever the oximeter is not in use.

### **Cleaning**

#### **Oximeter**

**WARNING:** Electrical shock and flammability hazard—Always turn the oximeter off and disconnect it from AC mains power before cleaning it.

#### **CAUTION:**

##### **Cleaning**

- Do not autoclave, pressure sterilize, or gas sterilize this oximeter.
  - Do not soak or immerse the monitor in any liquid.
  - Use the cleaning solution sparingly. Excessive solution can flow into the monitor and cause damage to internal components.
  - Do not touch, press, or rub the display panel with abrasive cleaning compounds, instruments, brushes, rough-surface materials, or bring it into contact with anything that could scratch the panel. Do not use solutions containing acetone or other harsh solvents to clean the display panel.
1. Clean the outer surface of the oximeter with a soft cloth dampened in a mild soap and water solution or isopropyl alcohol (70%).
  2. Allow the unit to dry completely before using it.

### Probes

**Important:** Follow the cleaning instructions on the user instructions shipped with the probe or sensor you are using.

In general, you may clean **reusable probes after each patient use**, as follows:

1. Disconnect the probe from the patient and the oximeter.
2. Clean with a soft cloth using a mild soap and water solution, or an isopropyl alcohol (70%) swab.
3. Allow the probe to dry completely before returning it to operation.

### Recharging the battery

**WARNING: Electrical shock hazard—Use only the charger supplied for this oximeter. Damage to the equipment may also result.**

**CAUTION:** Do **not** turn the oximeter on after the RECHARGE BATT alarm message has been displayed without first connecting it to the charger which is plugged into AC mains power. Operating the oximeter on battery power during a RECHARGE BATT alarm condition may permanently damage the lead-acid battery.

To recharge the battery, plug the charger into AC mains power and into the charging connector on the rear panel of the oximeter. Under normal conditions, the battery lasts for several hundred "charge-discharge" cycles. To obtain maximum battery life, recharge the battery pack whenever it is not in use. The battery will not overcharge.

**Note:** During the recharging process, the oximeter may be operated when it is plugged into AC mains power.

### Alarm messages

**LO BAT** appears when battery capacity is below 50%. Recharge the battery promptly when this message appears. **Note:** This message replaces pulse rate data readings on the display.

**RECHARGE BATT** appears for approximately 10 seconds, the audible alarm sounds, and the oximeter automatically shuts off. Immediately connect the charger into the oximeter and into AC mains power.

### Repair policy and procedure

Warranty repair and service must be performed by an Ohmeda Service Representative or at the Ohmeda Service and Distribution Center. To contact an Ohmeda Service Representative, call the nearest Ohmeda Regional Service Office listed on the back cover.

Do **not** use malfunctioning equipment. Make all necessary repairs or have the equipment repaired by an Ohmeda Service Representative. Parts listed in the service manual for this product may be repaired or replaced by a competent, trained person who has experience in repairing devices of this nature. After repair, test the equipment to verify that it complies with the published specifications.

**WARNING:** Do not remove the cover of the monitor. An operator may only perform maintenance procedures specifically described in this manual. Refer servicing to qualified service personnel trained in the repair of this equipment.

**CAUTION:** Detailed information for extensive repairs is included in the *Ohmeda 3760 Pulse Oximeter Service Manual* solely for the convenience of users who have proper knowledge, tools, and test equipment, and for service personnel trained by Ohmeda.

## **Obtaining service**

Contact the nearest Ohmeda Regional Service Office or Ohmeda representative listed on the back cover of this manual.

### **Packaging and return procedure**

If returning equipment to Ohmeda:

**Everyone:**

**Please clean contaminated/dirty equipment before returning.**

Package the equipment securely—in the original shipping container if possible—and enclose the following items:

1. A letter describing in detail any difficulties experienced with the Oximeter.
2. Warranty information--a copy of the invoice or other applicable documentation must be included.
3. Purchase order number to cover repair if the oximeter is not under warranty.
4. Ship-to and bill-to information.
5. Person (name, telephone/Telex/fax number, and country) to contact for any questions about the repairs.

When Ohmeda's warranty is not applicable, repairs are made at Ohmeda's current list price for replacement part(s) plus a reasonable labor charge.

#### **Inside the USA:**

**First**, call the Ohmeda Service and Distribution Center (800 999 6277) for instructions on your specific product and then ship it prepaid to the following address:

Ohmeda Service and Distribution Center  
7750 The Bluffs NW  
Austell, GA 30001  
Attention: Service Center

#### **Outside the USA:**

Send to your local authorized service office as shown on the back cover of this manual.

### Accessories

For a current list of Ohmeda probes, sensors, and their accessories, contact your Ohmeda representative.

Description	Part number
3760 Pulse Oximeter Operation and Maintenance Manual	6050-0002-840
3760 Pulse Oximeter Service Manual	6050-0002-853
Thermal printer paper	
Single:	0380-1500-034
Box of 5:	0380-1500-021

## **A/Specifications**

**Note:** All specifications are nominal and subject to change without notice.

### **SpO<sub>2</sub> accuracy (1 standard deviation)**

<b>SpO<sub>2</sub> range</b>	<b>Accuracy</b>	<b>Number of data points</b>
60% – 100%	2.4%	616
90% – 100%	1.5%	183
80% – 89.9%	2.1%	197
Below 60%	Unspecified	

These accuracy measurements are statistically derived and correlated to simultaneous oxygen averaging readings measured on an Ohmeda Blox 3700 Pulse Oximeter.

### **Interfering substances**

Carboxyhemoglobin may erroneously increase readings. The increase is approximately equal to the amount of carboxyhemoglobin present.

Dyes, or any substances containing dyes that change usual arterial pigmentation may cause erroneous readings.

See references at the end of this appendix.

### **Pulse rate accuracy (1 standard deviation)**

± 1.7% of current reading (assuming a constant pulse rate)

### **SpO<sub>2</sub> range**

0% to 100%

### **Pulse rate range**

40 to 235 beats per minute

### **Alarm limits**

Low SpO<sub>2</sub> = 80, 85, 90 (default), OFF

### **Oximeter**

Height:	8.1 cm	(3.2 in)
Width:	26.1 cm	(10.3 in)
Depth:	27.4 cm	(10.8 in)
Weight:	2.8 kg	(6.0 lb) without Printer
	3.2 kg	(7.0 lb) with Printer

### **Battery**

Sealed lead-acid, 12 volt, 1.9 Ampere-Hours

Charge Time (unit not operating, from a 25% charge capacity):

- 80% capacity in 4 hours
- 100% capacity in 4.5 hours

Operation time (from a fully charged battery to automatic shut-off at 25% capacity):

- Without Printer: approx. 7 hours
- With Printer: approx. 3.5 hours

### **Low battery indicator**

LO BAT message appears when battery capacity is below 50%

### **Electrical requirements**

Oximeter input power maximum: 15.5 volts AC  $\pm 5\%$ , 50 or 60 Hz (non-interchangeable) at 1.4 A

#### **100V charger**

Input (AC mains): 100 VAC  $\pm 10\%$  (217 mA max. current draw)

Output: 15.5 VAC at 1.4 A, 50/60 Hz

#### **120V charger**

Input (AC mains): 120 VAC  $\pm 10\%$  (180 mA max. current draw)

Output: 15.5 VAC at 1.4 A, 60 Hz

#### **220/240V charger**

Input (AC mains): 220/240 VAC  $\pm 10\%$  (100 mA max. current draw)

Output: 15.5 VAC at 1.4 A, 50 Hz

### **Environmental tolerances**

Operating temperature:	10° – 45° C	(50° – 115° F)
------------------------	-------------	----------------

Storage temperature (without paper):	-20° – 60° C	(-4° – 140° F)
--------------------------------------	--------------	----------------

Storage temperature (with paper):	-20° – 45° C	(-4° – 115° F)
-----------------------------------	--------------	----------------

Altitude:	to 3048 m	(10,000 ft)
-----------	-----------	-------------

Humidity:	5% to 80% non-condensing
-----------	--------------------------

**Note:** At temperature extremes, the LCD read-out may exhibit reduced contrast, ghosting or darkening. When returning from temperature extremes, allow the oximeter to return to room temperature before use.



**International Electrotechnical Commission classifications**

*Type of protection against electric shock:* Class I/Internal electrical power source.

*Degree of protection against electric shock:* Type BF

*Degree of protection against ingress of liquids:* Ordinary

*Mode of operation:* Continuous

*Recommended methods of sterilization or disinfection:* See 4/Maintenance and Service in this manual and appropriate section in the instructions shipped with your probes for recommended procedures for cleaning this equipment.

*Degree of safety of application in the presence of a flammable anesthetic mixed with air or with oxygen or nitrous oxide:* Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

**Printer**

Printer head life expectancy: approximately 500,000 lines

Paper (thermal-type) approximate dimensions:

Width:	5.7 cm	(2.25 in)
Roll diameter:	3.8 cm	(1.5 in)
Roll length:	15.8 cm	(52 ft)

One full roll of paper lasts approximately:

	60 Hz	50 Hz
Graph fast	22 hrs	17 hrs
Graph slow	67 hrs	62 hrs
Print fast	10 hrs	8 hrs
Print slow	45 hrs	44 hrs

## **References**

Barker, SJ, Ph.D., M.D.; Tremper, KK, Ph.D., M.D.; Hufstadler, S, M.D.; Hyatt, J; Zaccari, J. "The Effects of Carbon Monoxide Inhalation on Noninvasive Oxygen Monitoring" *Anesthesia Analg*, Vol 65:S12, 1986.

Barker, SJ, Ph.D., M.D.; Tremper, KK, Ph.D., M.D.; Hufstadler, S, M.D.; Hyatt, J; Zaccari, J. "Effects of Methemoglobinemia on Pulse Oximetry and Mixed Venous Oximetry" *Anesthesiology Supp.*, Vol. 67:3A, 1987.

Scheller, MS, M.D., Unger, RJ, M.D. "The Influence of Intravenously Administered Dyes on Pulse Oximetry Readings" *Anesthesiology*, Vol 65:A161, 1986.

Shippy, Merrill B., BS RRT; Petterson, Michael T., BA RRT CCPT; Whitman, Robert A., BA RRT CCPT; Shivers, Constance R., BA. "A Clinical Evaluation of the BTI Biox II Ear Oximeter" *Respiratory Care*, Vol 29:730, 1984.

Sidi, A, Rush, WR, Paulus, DA, M.D.; Gravenstein, N, M.D.; Davis, RF, M.D. "Effect of Fluorescein Indocyanine Green, and Methylene Blue on the Measurement of Oxygen Saturation by Pulse Oximetry" *Anesthesiology*, Vol 65:A132, 1986.

- accessories 4-4
- alarms
  - Alarm Limit key 1-4
  - Alarm silence key 1-4
  - alarm limits A-1
  - alarm silencing 1-4
  - low SpO<sub>2</sub> alarm limit violation 3-1
  - printer 2-4, 3-3
  - probe alarm conditions 3-1
  - signal and device failure problems 3-3
- battery
  - LO BAT message 3-3, 4-1, 4-2, A-2
  - RECHARGE BATT message 4-1, 3-3, 4-2
  - recharging 4-2
  - specifications A-2
- carboxyhemoglobin 1-1, 1-3, A-1
- cleaning
  - oximeter 4-1
  - probes 4-2
- clock set mode 2-3
- closing statistics 2-8, 2-9
- components, functional 1-3
- connectors
  - charger 1-5
  - probe plug 1-4
- data validity 2-2
- display LCD 1-4
- dyes (effect of) A-1
- features 1-4
- filled graph 2-7
- functional components 1-3
- graphic format 2-6
- interfering substances A-1
- International Electrotechnical Commission classifications A-3
- keys
  - Alarm Limit 1-4
  - Alarm silence 1-4
  - Power 1-4
  - for clock set mode 2-3
- maintenance
  - cleaning 4-1
  - recharging the battery 4-2
  - repair policy and procedure 4-2
  - obtaining service 4-3
- numeric format 2-5
- open graph 2-7
- oximeter
  - cleaning 4-1
  - electrical requirements A-2
  - features and controls 1-3
  - functional components 1-4
  - general operation guidelines 2-1
  - introduction 1-1
  - preoperative checklist 2-1
  - principles of operation 1-1
  - specifications A-2
- packaging and return procedure 4-3
- parts list, accessories 4-4
- precautions
  - cautions 1-7
  - warnings 1-5
- preoperative checklist 2-1
- visual inspection 2-1
- functional inspection 2-2
- principles of operation 1-1
- printer operation 2-4
  - alarms while printing 2-4, 3-3
  - closing statistics 2-8, 2-9
  - filled graph 2-7
  - graphic format 2-6
  - graphic resolution 2-6
  - interrupting printing 2-9
  - introduction 2-4
  - loading printer paper 2-10
  - numeric format 2-5
  - open graph 2-7
  - parts of a print-out 2-7
  - side panel—printer door 2-4
  - specifications A-3
  - stopping printing 2-9
  - switching formats 2-8
  - switching speeds 2-8
  - trend statistics 2-8, 2-10
  - thermal printer paper 2-4, 4-4
- probes
  - principles of operation 1-3
  - probe alarm conditions 3-1
  - preoperative checklist 2-1
- probe alarm messages
  - ID ERROR 3-2
  - LOW LIGHT 3-2
  - NO PROBE 3-1, 3-2
  - OFF PATIENT 2-2, 3-1, 3-2
  - PROBE/UNIT FAIL 3-2
- probe plug connector 1-4
- pulse rate specifications A-1
- repair policy and procedure 4-2
- references A-4
- service
  - obtaining service 4-3
  - repair policy and procedure 4-2
  - packaging and returning equipment 4-3
  - signal and data assurance -- SEE data validity 2-2
  - specifications

- alarm limits A-1
- battery A-2
- electrical requirements A-2
- environmental tolerances A-2
- interfering substances A-1
- International Electrotechnical Commission  
classifications A-3
- low battery requirements A-2
- 100V charger A-2
- 120V charger A-2
- oximeter A-1
- printer A-3
- pulse rate accuracy A-1
- pulse rate range A-1
- SpO<sub>2</sub> accuracy A-1
- SpO<sub>2</sub> range A-1
- 220/240V charger A-2
- statistics
  - closing 2-8, 2-9
  - trend 2-8, 2-10
- theory 1-1
- trend statistics 2-8, 2-10
- troubleshooting 3-3
- warranty W-1

## **Warranty**

The Ohmeda 3760 Pulse Oximeter is sold by Ohmeda only under the warranties set forth in the following paragraphs. Such warranties are extended only with respect to the purchase of this product directly from Ohmeda's Authorized Dealers as new merchandise and are extended to the first Buyer thereof, other than for resale.

### **Limited warranty**

Ohmeda warrants that the Ohmeda 3760 Pulse Oximeter meets its published specifications at the time of shipment from the factory. Ohmeda further warrants that this product was calibrated at the factory prior to shipment [using standard gravimetrically determined calibration gases. The calibration gases used to calibrate the instrument are warranted by Ohmeda.]

### **Products not under warranty**

The following items are not covered under this warranty: disposable items, service kits, and replacement parts. These items may be covered under a separate warranty. Consult Ohmeda for details.

### **Duration**

This product, other than its expendable parts, is warranted against functional defects in materials and workmanship for a period of twelve (12) months from the date of delivery to the user (in no event for a period of more than twelve (12) months from the date of original delivery by Ohmeda to an Ohmeda Authorized Dealer). If any part of this product proves defective under proper and normal use within the warranty period, and the proper planned maintenance procedures have been followed, as the purchaser's exclusive remedy, Ohmeda will repair or replace, at its sole discretion, the product or any defective part provided it is returned to Ohmeda Service within 30 days.

### **Limitation**

Ohmeda may at any time discharge its warranty obligation by repairing and returning the product to original factory performance. This may be accomplished by installing new or remanufactured assemblies or other repairs deemed appropriate by Ohmeda. Ohmeda's choice of repair or replacement shall be the sole remedy of the buyer or user.

### **Conditions**

This warranty is valid only when Ohmeda-trained personnel have performed installation and service on the product and all recommended planned maintenance procedures have been completed during the warranty period. Damage caused by the abuse or misuse of the product is not covered by this warranty. Ohmeda shall not be liable for damage resulting from the improper installation or the misuse of the product.

### **Exclusion of warranties**

Oral statements about the product do not constitute warranties, shall not be relied on by the buyer or user, and are not part of any warranty extended by Ohmeda.

**Except as set forth in this limited warranty, Ohmeda makes no warranties, expressed or implied, including the implied warranty of merchantability and the implied warranty of fitness for a particular purpose. Except for the obligations under this limited warranty, Ohmeda shall not have any obligation or liability for any incidental or consequential damages (including those from commercial loss) or other loss, damage, or injury resulting directly or indirectly for the product.**

## North America

### United States

#### **USA Customer Service**

Ohmeda Inc  
Ohmeda Drive  
PO Box 7550  
Madison WI 53707 7550  
USA  
Tel 800 345 2700  
Fax 608 221 4384

#### **USA Technical Support**

Ohmeda Inc  
Ohmeda Drive  
PO Box 7550  
Madison WI 53707 7550  
USA  
Tel 800 345 2755  
Fax 608 221 3618

#### **USA Service and Distribution Center**

Ohmeda Inc  
7750 The Bluffs NW  
Austell GA 30001  
USA  
Tel 800 345 2755  
Tel 770 739 4774  
Fax 770 739 4770

### Canada

#### **Customer Service and Service Centre**

Ohmeda  
5975 Falbourne St., Unit 2  
Mississauga Ont  
Canada L5R 3V8  
Tel 905 568 9533  
Fax 905 568 9799  
Telex 06989362

### Asia/Pacific

#### India

Ohmeda  
B.O.C. India Ltd  
Oxygen House  
P43 Taratala Road  
Calcutta 700 088  
India  
Tel 91 33 478 4342  
Fax 91 33 478 4410

## Japan

Ohmeda Ltd  
Technical Center  
No. 21 Maruya YA Bldg.  
2 16 11 Kaname-cho  
Toshima-ku  
Tokyo 171 Japan  
Tel 81 3 3973 4055  
Fax 81 3 3973 9102

### Malaysia

Ohmeda  
Malaysian Oxygen Berhad  
13 Jalan 222  
46100 Petaling Jaya  
Selangor  
Malaysia  
Tel 60 3 755 4233  
Fax 60 3 755 6389

### Singapore

Ohmeda Pte Ltd  
No 198 Yishun Avenue 7  
Singapore 768926  
Tel 65 753 6311  
Fax 65 753 1817

### Australia

Ohmeda Sales and Service  
1/2 149 Arthur Street  
PO Box 356  
Homebush NSW 2140  
Australia  
Tel 61 132 229  
Fax 61 297 461 796

## Europe

### Austria

Ohmeda Austria GmbH  
Am Concorde Park E1/7  
A-2320 Schwechat  
Austria  
Tel 43 1 706 36 60  
Fax 43 1 706 36 60 11

### Belgium

Ohmeda  
Divisie van BOC NV  
Maaltecenter Blok G  
Derbystraat 251  
9051 Gent  
Belgium  
Tel 32 9 220 4757  
Fax 32 9 221 8541

## Denmark

BOC Ohmeda A/S  
Carl Gustavs Gade 3  
Høje-Taastrup  
DK-2630 Taastrup  
Denmark  
Tel 45 43 71 43 72  
Fax 45 43 71 72 71

### Finland

BOC Ohmeda Oy  
PL 78 Kuunkehrrä 2 B  
02211 Espoo  
Finland  
Tel 358 9 803 73 44  
Fax 358 9 803 05 69

### France

Ohmeda  
Parc de Pissaloup  
8 Avenue Jean d'Alembert  
78190 Trappes-Cedex  
France  
Tel 33 1 30 68 60 00  
Fax 33 1 30 68 60 01

### Germany

Ohmeda  
Langemarckplatz 3  
Postfach 2880  
D-91054 Erlangen  
Germany  
Tel 49 91 31 81 01 10  
Fax 49 91 31 81 01 77

### Italy

Ohmeda  
Via Carpaccio 33  
20090 Trezzano sul  
Naviglio  
Milan  
Italy  
Tel 39 2 445761  
Fax 39 2 48403233  
Telex 314479 I

### Netherlands

Ohmeda BV  
Rembrandtlaan 1A  
3723 BG Bilthoven  
Netherlands  
Tel 31 30 2289 289  
Fax 31 30 2288 366

## Spain

Ohmeda  
C/Manuel Tovar 26  
28034 Madrid  
Spain  
Tel 34 1 334 26 00  
Fax 34 1 358 12 84

### Sweden

BOC Ohmeda AB Nordic  
PO BOX 631  
Östra Sandgatan 12  
S 251 06 Helsingborg  
Sweden  
Tel 46 42 17 87 40  
Fax 46 42 17 87 47

### United Kingdom

Ohmeda  
Ohmeda House  
71 Great North Road  
Hatfield Hertfordshire  
AL9 5EN England  
Tel 44 1707 263570  
Fax 44 1707 260191  
Telex 915128

### Latin America Caribbean

Ohmeda  
9200 S Dadeland Blvd  
Suite 610  
Miami FL 33156  
USA  
Tel 305 670 9938  
Fax 305 670 4834

### Middle East

#### Dubai

Ohmeda  
Middle East Operations  
PO Box 5527  
Dubai  
United Arab Emirates  
Tel 97 14 822653  
Fax 97 14 822659