

Power Output for ArthroCare System Quantum Controller

Ablation- Load - **250** ohms / unit of measure in Watts

Setting	Acceptable Range (+/- 20%)	Typical	Actual
*****	*****	*****	*****
1	28 - 42	35	_____
2	45 - 67	56	_____
3	66 - 98	82	_____
4	90 - 134	112	_____
5	119 - 179	149	_____
6	152 - 228	190	_____
7	188 - 282	235	_____
8	230 - 344	287	_____
9	274 - 400 (Max)	343	_____

Serial # _____

Control # _____

Date _____

Coag - Load - **250** ohms / unit of measure in Watts

Setting	Acceptable Range (+/- 20%)	Typical	Actual
*****	*****	*****	*****
1	12 - 18	15	_____
2	27 - 41	34	_____

Ablation- Load - **300** ohms / unit of measure in Watts

Setting	Acceptable Range (+/- 20%)	Typical	Actual
*****	*****	*****	*****
1	23 - 35	29	_____
2	37 - 55	46	_____
3	54 - 82	68	_____
4	76 - 114	95	_____
5	100 - 150	125	_____
6	128 - 192	160	_____
7	158 - 236	197	_____
8	192 - 288	240	_____
9	230 - 344	287	_____

Coag - Load - **300** ohms / unit of measure in Watts

Setting	Acceptable Range (+/- 20%)	Typical	Actual
*****	*****	*****	*****
1	10 - 14	12	_____
2	23 - 35	29	_____